

# DSLSA9 SIT-TO-STAND BUTTOCK STRAP: APPLICATION INSTRUCTIONS

#### PLEASE SAVE THESE INSTRUCTIONS FOR FUTURE USE

Note: The most current version of these instructions can be found online at www.grahamfield.com

- $\triangle$  WARNING: If components are damaged or missing, contact your dealer immediately.

## **APPLICATION**

- 1. Check the patient's weight. Ensure that the patient's weight does not exceed the Buttock Strap's maximum weight capacity.
- 2. Review the following procedures carefully before attempting lift; call Graham-Field Technical Support at 800-347-5678 with any questions.
- 3. Use long loops, center loops, or short loops as shown in Figure 2, to position Buttock Strap to provide the greatest support to the patient's middle and lower body during the following steps.


Figure 1


A

B

a

b







- 4. Attach Buttock Strap strap **A** to Sit-to-Stand Lift mast lower loop **a** as shown in Figures 1 and 3.
- 5. Bring Buttock Strap around the back of the patient's buttocks and attach Buttock Strap strap **B** to Sit-to-Stand Lift mast lower loop **b** as shown in Figures 1 and 3.
- 6. Do not lock caster brakes, in order to allow the lift to properly balance weight.
- 7. Proceed with lift and transfer as shown in Figure 4.

DSLSA9-INS-LAB-RevA09

#### **MAINTENANCE**

To clean Buttock Strap, use a mild soap and water solution and a clean cloth. Rinse with clear water, then allow to air dry.

▲ Caution: Never use solvents, abrasive cleaners, or scouring pads on any part of the Buttock Strap.

#### **WARRANTY**

GF Health Products, Inc. offers a one-year limited warranty against manufacturer's defects on the Sit-to-Stand Buttock Strap. If a product is deemed to be under warranty, GF Health Products, Inc. shall provide, at its option, (1) replacement of any defective part or product or (2) a credit of the original selling price made to GF Health Products, Inc.'s initial customer. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to GF Health Products, Inc.


DSLSA9

34 1/4"

### **FEATURES**

**Sit-to-Stand Buttock Strap** Designed for use with

LF2020 STS Lift. 400 lb maximum weight capacity. **DSLSA9** 1 ea



71 1/4"

5 1/2"

7 1/2"



GF Health Products, Inc. is not responsible for typographical errors. Packaging, warranties, products and specifications are subject to change without notice. Graham-Field and Lumex are registered trademarks of GF Health Products, Inc.

© October 2009 GF Health Products, Inc. • 2935 Northeast Parkway • Atlanta, GA 30360 • tel 800-347-5678 • fax 800-721-0601 • www.grahamfield.com

DSLSA9-INS-LAB-RevA09 2