Electrical Neuromuscular Stimulation (EMS), Digital Model GF-TX5EMS

Operation Manual

Read this manual before operating your GF-TX5EMS.

Save this manual for future use.

The most current version of this manual can be found online at
www.grahamfield.com.

CONTENTS

GENERAL DESCRIPTION	3
WHAT IS EMS?	3
INDICATIONS AND CONTRAINDICATIONS	3
SAFETY	4
ABOUT THE DEVICE.....	7
EXPLANATION OF KEY / KNOB CONTROL FUNCTIONS	8
ATTACHING THE LEAD WIRES.....	9
ELECTRODE SELECTION AND CARE.....	9
TIPS FOR SKIN CARE	9
CONNECTING THE DEVICE.....	10
BATTERY INFORMATION.....	11
CARING FOR YOUR DEVICE.....	12
TROUBLESHOOTING.....	12
TECHNICAL SPECIFICATIONS.....	13
LIMITED WARRANTY	15
MANUAL DE OPERACIÓN EN ESPAÑOL.....	17

GF Health Products, Inc. is not responsible for typographical errors. For the most updated and current information on packaging, warranties, products and specifications, including the most current version of these instructions, please visit our website at www.grahamfield.com.

Graham-Field and **Grafo** are registered trademarks of GF Health Products, Inc.

GENERAL DESCRIPTION

EMS, Electrical Neuromuscular Stimulation (Electrical Myostimulation), is a battery-operated pulse generator that sends electrical impulses to electrodes attached to the body to stimulate motor nerves and cause contraction and relaxation of muscles. It has proven valuable as a method of pain therapy and is of great assistance to the experienced therapist.

With some indications, physicians / clinicians can prescribe EMS units to patients for their use at home.

This unit is a dual-channel digital stimulator for active treatment application, which has a Liquid Crystal Display indicating operation modes and output, as well as an 8-bit microcomputer for controlling the system.

The electronics of the unit create electric impulses; the intensity, duration, frequency per second and modulation of these impulses can be adjusted.

WHAT IS EMS?

An injured muscle usually experiences little if any movement. EMS therapy remedies this by using low voltages to stimulate motor nerves to cause involuntary muscular contractions.

Like exercise, EMS helps strengthen the injured area, and has been found to effectively treat a variety of musculoskeletal and vascular conditions. Common candidates for EMS therapy are patients recovering from orthopedic surgery, muscle strains or tears, or athletes who have undergone cartilage or tendon repair. EMS is non-invasive and does not use pharmaceuticals.

INDICATIONS AND CONTRAINDICATIONS

Read the operation manual before using this EMS device.

Federal law (USA) restricts this device to sale by or on the order of a physician.

Observe your physician's / clinician's precise instructions and allow them to show you where to apply the electrodes. For a successful therapy, the correct application of the electrodes is an important factor. Carefully write down the settings your physician / clinician recommends.

Indications

- Relaxation of muscle spasm
- Prevention or retardation of disuse atrophy
- Increasing local blood circulation
- Muscle re-education
- Immediate postsurgical stimulation of calf muscles to prevent venous thrombosis
- Maintaining or increasing range of motion

Powered muscle stimulators should only be used under medical supervision for adjunctive therapy for the treatment of medical diseases and conditions.

Contraindications

 **WARNING: This device should not be used by patients with cardiac demand pacemakers.**

SAFETY

Always follow basic safety precautions, including the following:

 **WARNING: Indicates a potential hazard situation or unsafe practice that, if not avoided, could result in death or serious personal injury.**

 **CAUTION: Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor or moderate personal injury.**

 **NOTICE: Indicates a potential hazard or unsafe practice that, if not avoided, could result in product / property damage.**

Warnings

 **WARNING: This device does not have AAP/APG protection.**

Explosion hazard is possible if used in the presence of explosives, flammable materials or flammable anesthetics.

 **WARNING: The long-term effects of chronic electrical stimulation are unknown.**

-  **WARNING: Do not place electrodes over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.**
-  **WARNING: Do not place electrodes over the neck or mouth. This may result in spasms of the laryngeal and pharyngeal muscles, and the contraction may be strong enough to close the airway or cause difficulty in breathing.**
-  **WARNING: Do not apply stimulation transthoracically, because the introduction of electrical current into the heart may cause cardiac arrhythmia.**
-  **WARNING: Do not apply stimulation transcranially.**
-  **WARNING: Do not apply stimulation over swollen, infected or inflamed areas or skin eruptions; e.g. phlebitis, thrombophlebitis, varicose veins, etc.**
-  **WARNING: Do not apply stimulation over or in proximity to cancerous lesions.**
-  **WARNING: The safety of the device during pregnancy or delivery has not been established.**
-  **WARNING: Caution should be used for patients with suspected or diagnosed heart problems.**
-  **WARNING: Caution should be used when applying this device to patients suspected of having heart disease. Further clinical data is needed to show if there are adverse side effects on individuals with heart disease.**
-  **WARNING: Caution should be used for patients with suspected or diagnosed epilepsy.**
-  **WARNING: Caution should be used in the presence of the following:
When there is a tendency to hemorrhage following acute trauma or fracture;
following recent surgical procedures when muscle contraction may disrupt the healing process; over the menstruating or pregnant uterus; and over areas of the skin which lack normal sensation.**
-  **WARNING: Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.**

-  **WARNING: Keep this device out of the reach of children.**
-  **WARNING: Use this device only with the leads and electrodes recommended for use by GF Health Products, Inc.**
-  **WARNING: Do not use this device while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.**
-  **WARNING: EMS devices should be used only under the continued supervision of a physician / clinician.**
-  **WARNING: Notice for California Customers- California Proposition 65**
WARNING: This product contains a chemical known to the State of California to cause cancer and reproductive or developmental harm.

Precautions / Adverse Reactions

-  **CAUTION: Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. Using an alternate conductive medium or alternate electrode placement can usually reduce the irritation. Consult your physician / clinician before using an alternative conductive medium or electrode placement.**
-  **CAUTION: Isolated cases of skin irritation may occur at the site of electrode placement following long-term application.**
-  **CAUTION: If skin irritation occurs EMS treatment should be stopped and electrodes removed until the cause of the irritation can be determined.**
-  **CAUTION: Effectiveness is highly dependent upon patient selection of a doctor qualified in the management of pain patients.**
-  **CAUTION: If the device treatment becomes ineffective or unpleasant, stimulation should be discontinued until reevaluation by a physician / clinician.**
-  **CAUTION: Always turn the device OFF before applying or removing electrodes.**

ABOUT THE DEVICE

This device is a battery-operated device that includes two controllable output channels. This device creates electrical impulses whose intensity, duration, and modulation can be altered. The device controls are easy to use and the slide cover protects accidental changes in settings.


System Components

Your device will include the following components or accessories:

- EMS unit
- Carrying case
- Lead wires
- Operation manual
- Electrodes

Also required (not included): One 9-Volt alkaline or nickel-cadmium rechargeable battery.

Device Controls

Slide Cover

This cover located on the front of the unit conceals the controls for DOWN, MODE, SET, and UP. Press the top portion of the cover and pull down in order to open the cover.

EXPLANATION OF KEY / KNOB CONTROL FUNCTIONS

DOWN Key	This key decreases the setting of treatment parameters. *	This key regulates the number of pulse width, pulse rate, and cycle on/ramp/off time.
	*Decrease by pressing the DOWN key: The width can be adjusted by 10μ/step. The rate can be adjusted in 1Hz~20Hz by 1Hz/step, 20Hz~150Hz by 5Hz/step. The cycle on/ramp/off time can be adjusted by 1 second/step.	
MODE Key	This key selects <u>Timer Mode</u> , <u>Stimulation Mode</u> , <u>Pulse Width Mode</u> , <u>Pulse Rate Mode</u> , or <u>Cycle On/Ramp/Off Time</u>	This key changes the treatment parameter. Each time the MODE key is pressed, the next parameter will display. The selected treatment parameter in the current mode will flash.
SET Key	This key switches between the different settings within <u>Timer Mode</u> and <u>Stimulation Mode</u> .	Each time the SET key is pressed, the parameter will change to the next setting within the parameter. The selected setting will flash. When the desired setting is flashing, press the MODE key to confirm the choice.
UP Key	This key increases the setting of treatment parameters. *	This key regulates the number of pulse width, pulse rate, and cycle on/ramp/off time.
	*Increase by pressing the UP key: The width can be adjusted by 10μs/step. The rate can be adjusted in 1Hz~20Hz by 1Hz/step, 20Hz~150Hz by 5Hz/step. The cycle on/ramp/off time can be adjusted by 1 second/step.	
Ch1 / Ch2 Knobs	Channel 1 and Channel 2 Intensity Control Knobs	These knobs control the strength of the stimulation and also function as ON / OFF controls.

ATTACHING THE LEAD WIRES

 **WARNING:** Ensure the device is **OFF** before connecting the lead wires.

 **WARNING:** Never insert the lead wire plug into an AC power supply socket. Personal injury and damage to the EMS unit could occur.

 **NOTICE:** Use care when you plug and unplug the wires. Pulling on the lead wire instead of its insulated connector may cause wire breakage.

The lead wires provided with the device insert into the ports located on the top of the unit. If only one lead will be used, plug it into the channel 1 port. After connecting the wires to the unit, attach each wire to an electrode.

Lead wires provided with the device are compliant with mandatory compliance standards set forth by the FDA.

ELECTRODE SELECTION AND CARE

Using Electrodes

Use the electrodes as prescribed. Follow application procedures outlined in electrode packaging to maintain stimulation and prevent skin irritation. The electrode packaging provides instructions for care, maintenance, and proper storage of electrodes.

TIPS FOR SKIN CARE

Good skin preparation is important for effective and comfortable use of your EMS device.

- Always clean the electrode site with mild soap and water solution, rinse well, and dry thoroughly prior to any electrode application.
- Any excess hair should be clipped, not shaved, to ensure good electrode contact with the skin.
- If a skin treatment or preparation is recommended by your physician / clinician, apply the skin treatment as recommended, let dry, and apply electrodes as directed. Following these recommendations will both reduce the chance of skin irritation and extend the life of your electrodes.

- Avoid excessive stretching of the skin when applying electrodes. Proper application is best accomplished by applying the electrode, then smoothly pressing it in place from the center outward.
- When removing electrodes, always remove by pulling in the direction of hair growth.
- It may be helpful to rub skin lotion on electrode placement area when not wearing electrodes.

CONNECTING THE DEVICE

Inserting the Battery

Turn the device to the OFF position before inserting or removing the battery. When inserting the battery, ensure the battery polarity (+ and –) markings match the markings on the device.

Preparing the Skin

Prepare the skin as previously described and according to the instructions provided with your electrodes. Before attaching the electrodes, identify the area that your physician / clinician has recommended for electrode placement.

1. **Connect the lead wires to the electrodes:** connect the lead wires to the electrodes before applying the electrodes to the skin.

CAUTION: Ensure both Intensity Controls for Channel 1 and 2 are turned to the OFF position (counterclockwise) before applying the electrodes.

2. **Place electrodes on the skin:** place the electrodes on the skin as recommended by your physician / clinician.
3. **Insert lead wire connector into the device:** plug end of lead wire into the channel output port (jack) to be used; push the plug in as far as it will go.
4. **Select treatment settings:** ensure your unit is still set to the proper settings recommended by your physician / clinician.
5. **Adjusting Channel Intensity Control:** locate the Intensity Control Knob (Channel 1 or 2) at the top of the unit. Slowly turn the Intensity Control Knob clockwise until the stimulation is at the level recommended by your

physician / clinician. (If you don't feel anything, turn the Knob OFF then ON again and carefully turn the Control Knob until you feel a tingling or slight twitch under or around the electrodes.) Always start with the lowest setting and increase the intensity slowly.

If the stimulation levels are uncomfortable or become uncomfortable, reduce the stimulation intensity to a comfortable level; or cease stimulation and contact your physician.

6. Setting the Patient Compliance Counter:

- a) To turn on the Patient Compliance Counter: While the unit is ON, hold down the UP key and press the MODE key simultaneously.
- b) To reset the Patient Compliance Counter: While the unit is ON, press the UP key and the MODE key simultaneously (this will take you into the Compliance Counter), then push the DOWN key and press the MODE key simultaneously.
- c) Press the UP and MODE keys simultaneously to return to the treatment status.

BATTERY INFORMATION

When the LCD (Liquid Crystal Display) low battery mark illuminates, the battery has become too weak to power the unit, and the existing battery should be replaced with a new battery. At this point, the unit will turn OFF until a new battery is inserted.

▲ NOTICE: GF Health Products, Inc. recommends the use of only a 9V alkaline or nickel-cadmium rechargeable battery with this device.

Replacing the Battery

When the LCD low battery mark illuminates, and the unit does not remain illuminated once turned on, the battery should be replaced.

1. Turn unit OFF.
2. Remove the front panel cover by pressing on the top of the panel and pressing down in order to slide the panel down. Continue sliding the panel downwards until the panel is completely removed from the unit. This will reveal the battery compartment.

3. Remove the discharged battery from the device.
4. Place new battery in the compartment. Note: Be sure the proper polarity (+ and –) markings match the markings in the device.
5. Replace the front panel cover.
6. Dispose of the old battery according to local guidelines and regulations.

CARING FOR YOUR DEVICE

Your device may be cleaned by wiping gently with a damp cloth moistened with mild soap and water. Do not immerse the device in water or other liquids.

Wipe lead wires with damp cloth moistened with soap and water. Do not immerse the lead wires.

To properly store the device for an extended period of time, remove the battery from the unit. Place the unit and accessories in the carrying case provided and store in a cool, dry location.

TROUBLESHOOTING

If the device does not function properly:

1. Ensure the battery is properly installed or replace the battery. Be sure to observe proper polarity markings when replacing the battery. If the LCD low battery mark illuminates when the unit is turned on, replace the battery and check again.
2. If the intensity has been adjusted and no stimulation is felt, check to ensure the lead wires are properly connected and the electrodes are properly applied to the skin. If the unit appears to be functioning and no stimulation is felt, the lead wires or electrodes may need to be replaced.
3. If the battery appears to be charged and the unit is not functioning, turn both Intensity Control Knobs to the OFF position (counterclockwise). Then gradually turn the Intensity Control Knob (clockwise) until stimulation is felt. If device still is not working, turn the unit off and contact your authorized GF Health Products, Inc. distributor.


If there is any other problem, please contact an authorized GF Health Products, Inc. distributor. Do not try to repair a defective device.

TECHNICAL SPECIFICATIONS

Channel	Dual, isolated between channels	
Pulse Intensity	Adjustable 0 mA – 100 mA peak into 500 ohm load each channel, constant current	
Pulse Width	30µs~260µs, adjustable by 10µs/step	
Pulse Rate	1Hz~150Hz, adjustable	1Hz~20Hz adjustable by 1Hz/step
		20Hz~150Hz adjustable by 5Hz/step
Timer	10, 20, 30, 40, 50, 60 minutes and continuous	
Patient Compliance Counter	To show the usage time by user, counted per each 30 minutes treatment time	
Patient Lock System	To prevent changes in any set parameter by physician / clinician. To unlock or lock the unit, press the UP and DOWN keys simultaneously for 2-3 seconds	
Wave Form	Symmetrical bi-phasic square pulse	
LCD	Displays Timer, Function Mode, Pulse Width, Pulse Rate, Cycle ON/ Ramp/OFF Time, Low Battery, Patient Compliance Counter, and Patient Lock System	
Function Mode	C	Constant Mode
	S	Synchronous Mode
	A	Alternating Mode
	Cycle ON Time (Contraction)	Adjustable 1-30 second, 1 second/step
	ON Ramp	Adjustable 1-10 second, 1 second/step
	Cycle OFF Time (Relaxation)	Adjustable 1-30 second, 1 second/step
Voltage	50 volt (peak to peak), at 500 ohm load	
Power Source	9 Volt battery	
Battery Life	Approximately 70 hours at normal setting	

Note: All electrical specifications are $\pm 10\%$ 500 ohm load.

LCD (Liquid Crystal Display)

The LCD displays timer (10, 20, 30, 40, 50, 60 minutes and C continuous), function modes (S, C, A, Cycle ON / Ramp / OFF time), pulse width, pulse rate, low battery, and patient compliance counter.

1. To check the LCD function, turn device ON and all the parameters will be displayed for 1~2 seconds.
2. After 1~2 seconds, the LCD will go to previous operation parameter.

Changing the LCD Parameters

- To set the treatment parameters on the EMS device, the unit must be in the UNLOCK position. If the device is in the LOCK position it, can be unlocked by pressing the UP and DOWN keys simultaneously for 2~3 seconds.
- To switch between the different treatment parameters, press the MODE key.
- Once in the desired mode, press the UP, DOWN, or SET key until the desired treatment parameter is obtained.
- Once all of the desired treatment settings are displayed, the unit may be locked by pressing the UP and DOWN keys simultaneously for 2-3 seconds while in the timer function.

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a component is warranted by a third party, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted.

This limited warranty shall only apply to defects that are reported to GF's customer service team within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. This limited warranty is not transferable.

The warranted components and time period are set forth below:

GF-TX5EMS Electrical Neuromuscular Stimulation (EMS): one year

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

This limited warranty shall only apply to defects that are reported to the Distributor from whom the Customer purchased the product within the applicable warranty period. If there is not a Distributor, you must contact GF directly by calling 770-368-4700, sending a fax request to 770-368-2386, or by e-mailing a request to cs@grahamfield.com. Specific directions will be provided by the Customer Service Representative. Failure to abide by the specific directions will result in denial of the warranty claim.

EXCLUSIONS

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a consumable nature including, but not limited to: casters, filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 6) Any labor or shipping charges incurred in the replacement part installation or repair;
- 7) Costs and expenses of regular maintenance and cleaning; and
- 8) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS.

THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

NOTES:

- 1) Additional terms and conditions may apply.
- 2) Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within thirty (30) days of the invoice date.

Manufactured for:
GF Health Products, Inc.
2935 Northeast Parkway
Atlanta, Georgia 30360 USA

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© 2010 GF Health Products, Inc.