



GF6650A-1: HOMECARE CLAMP-ON HALF RAIL INSTALLATION & OPERATION INSTRUCTIONS

SAVE THESE INSTRUCTIONS FOR FUTURE USE

Info: The most current version of these instructions can be found online at:

www.grahamfield.com

- ⚠ WARNING: Important! Read and understand these instructions before installing or using the GF6650A-1 Homecare Clamp-On Half Rail. If you do not understand any part of these warnings, cautions or instructions, contact a healthcare professional for direction in the use of this product. If the Homecare Clamp-On Half Rail is not properly installed and adjusted, personal injury and damage to the half rail could result.**


For information on Bed Rail safety, see the following brochure: "A Guide to Bed Safety: Bed Rails in Hospitals, Nursing Homes and Home Health Care: The Facts" on the FDA website, at:

www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/GeneralHospitalDevicesandSupplies/HospitalBeds/

- ⚠ WARNING: The Homecare Clamp-On Half Rail is designed for use in the homecare setting. The Homecare Clamp-On Half Rail is not intended for institutional use.**
- ⚠ WARNING: If components are damaged or missing, contact your dealer immediately. DO NOT use substitute parts. Use only Lumex replacement parts. Non-Lumex replacement parts could cause entrapment, personal injury, and damage to the half rail.**
- ⚠ WARNING: The Homecare Clamp-On Half Rail is intended to define the edge of the sleep surface and to help prevent an individual from inadvertently rolling out of bed. DO NOT use the half rail as a restraint.**
- ⚠ WARNING: Although it has no specific weight limitation, the Homecare Clamp-On Half Rail may deform or break when subjected to excessive side pressure.
DO NOT exert side pressure on the half rail.
DO NOT use it as a push handle for moving the bed.
DO NOT use it as an aid in transfer.
These activities could result in personal injury and damage to the half rail.**
- ⚠ WARNING: Notice for California Customers - California Proposition 65 WARNING: This product can expose you to chemicals including DEHP which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov/product.**
- ⚠ WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper installation or use of this product.**

INSTALLATION

Movable clamp disassembly




The Homecare Clamp-On Half Rail, as shipped, is shown at above left. Before installation, you must disassemble the movable clamp.


Unscrew the knob from the movable clamp assembly to separate the knob, spacer, and threaded clamp, as shown at above left.

Homecare Clamp-On Half Rail installation

Info: The Homecare Clamp-On Half Rail may be installed on either head deck angle or foot deck angle. The Homecare Clamp-On Half Rail CAN NOT be installed on pan deck.




install movable clamp jaw
over deck angle




install fixed clamp jaw
under deck angle

1. Attach movable clamp jaw over deck angle, with threaded rod facing outward, as shown at above left.
2. Position half rail with fixed clamp jaw toward bed and rails facing outward, as shown at above right. Position half rail mounting slot over movable clamp threaded rod and slide until fixed clamp jaw is installed under deck angle.



replace spacer on movable clamp



replace knob on movable clamp


3. Install spacer over end of movable clamp threaded rod and slide it against half rail flat surface, as shown at above left.
4. Install knob on end of movable clamp threaded rod, as shown at above right. Screw on until the clamp that secures the half rail to the bed frame is tight.

RAISE AND LOWER HOMECARE CLAMP-ON HALF RAIL

See Homecare Clamp-On Half Rail picture at right.

To raise half rail: Pull out retainer bar, raise rail to highest position, and release retainer bar, ensuring that rail is securely trapped by retainer hooks.

To lower half rail: Raise rail slightly to release retainer hooks, pull out retainer bar, lower rail, and release retainer bar.



raise and lower half rail

- ⚠ **WARNING:** Before use, ensure the Homecare Clamp-On Half Rail is assembled and installed as described, all attaching hardware is securely tightened, and the Homecare Clamp-On Half Rail is securely in place.
- ⚠ **WARNING:** During use, inspect frequently to ensure all attaching hardware is securely tightened and the Homecare Clamp-On Half Rail is securely in place.

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a component is warranted by a third party, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted. This limited warranty shall only apply to defects that are reported to GF's customer service team within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. This limited warranty is not transferable. Within the guidelines set forth in this document, this product is warranted for three (3) years. The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

This limited warranty shall only apply to defects that are reported to the Distributor from whom the Customer purchased the product within the applicable warranty period. If there is not a Distributor, you must contact GF directly by calling 1-770-368-4700, sending a fax request to 1-770-368-2386, or by e-mailing a request to cs@grahamfield.com. Specific directions will be provided by the Customer Service Representative. Failure to abide by the specific directions will result in denial of the warranty claim.

EXCLUSIONS

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a non-durable nature including, but not limited to: casters, filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 6) Any labor or shipping charges incurred in the replacement part installation or repair;
- 7) Costs and expenses of regular maintenance and cleaning; and
- 8) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

- 1) Additional terms and conditions may apply.
- 2) Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within three (3) days of the invoice date.



1.770.368.4700

Information contained herein is subject to change.
The most current and complete product information can be found on our website.
www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex and Graham-Field are trademarks of
GF Health Products, Inc., One Graham-Field Way, Atlanta GA 30340-3140
GF Health Products, Inc. is an ISO 13485 Certified Company.
Manufactured for GF Health Products, Inc.
Made in China