

Oxygen Regulators

- High-pressure regulators with flow control valves, designed to maintain a pressure of 50 psi and provide a variable flow rate within 12 settings when provided with a gaseous oxygen supply from 500 to 2200 psi
- New click style flow control
- Bumper-protected contents gauge
- Five Year Warranty
- Designed to meet the requirements of HCPCS Code: E1353

Adult

JB0150-081	0-8 LPM Yoke style/CGA870	1 ea
JB0150-082	0-15 LPM Yoke style CGA870	1 ea
JB0150-083	0-8 LPM Yoke style/CGA870 (DISS fitting)	1 ea
JB0150-085	0-8 LPM H-Nut style/CGA540	1 ea
JB0150-086	0-15 LPM H-Nut style/CGA540	1 ea



H-Nut style



Yoke style