

**AltaDyne<sup>™</sup> Alternating Pressure/Low Air Loss Mattress**

- An economical, portable air flotation system offers both alternating pressure and low air loss
- Intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort
- Group II mattress system suitable for use in the prevention and treatment of pressure ulcers
- Non-alternating, zoned (sacral region) low air loss can be selected by simply pressing the “static” switch
- The pressure control knob allows the firmness of the mattress to be adjusted for patient comfort
- LED indicator lights for low pressure alert
- Lighted ON/OFF indication switch
- The flotation mattress consists of individual transverse air cells designed for maximum support and pressure distribution
- Quilted nylon PU cover
- Quick CPR deflation is accomplished by simply disconnecting the quick-connect hose couplings from the control unit
- Electrical rating: 120 VAC/ 60Hz
- UL60601-1 and CSA C22.2 listed
- Air input > 4 lpm
- Cycle time: 8 minutes
- 350 lb maximum weight capacity, EVENLY DISTRIBUTED
- Pocket for optional 2” foam base
- Designed to meet the requirements of HCPCS code: E0277
- 1 year limited warranty
- Latex Free



<b>750000</b>	78.7" x 35.4" x 8"	1 ea
<b>763400F</b>	Optional 2" Foam Base	1 ea

