

Fast Alert Basic Patient Alarm

- Alarm notifies staff or caregiver when user gets out of bed or chair
- Pressure sensitive pad with 1 year warranty connects to alarm
- Alarm easily secures with self contained clip onto areas such as wheelchairs
- Convenient dual volume settings: 70dB and 90dB
- Delay function allows for 0 or 2 second delay to reduce likelihood of false alarms with restless users
- Requires one 9V battery (not included)
- Comes with protective boot cover to prevent damage if dropped
- Off switch and Low battery indicator
- 1 Year Limited Warranty
- Pad size:
 Chair - 10" x 15"
 Bed - 10" x 29"

GF13701	Basic Patient Alarm only	1 ea
GF13701B	Basic Patient Alarm with Bed Pad	1 set
GF13701C	Basic Patient Alarm with Chair Pad	1 set

Replacement Parts

106375	Pressure Sensor Bed Pad (1 Year Warranty)	1 ea
106380	Pressure Sensor Chair Pad (1 Year Warranty)	1 ea
GF1370-BCLIP	Bed Clip	1 ea
GF1370-COVER	Replacement Protective Boot Cover	1 ea

