

Versa-Helper Trapeze

- Helps patient safely change position in bed and transfer from bed to other equipment with minimal assistance.
- Rubber bumper guards protect walls.
- Handbar location fully adjustable.
- Height adjustable from 43" to 52".
- VersaGuard-coated clamps and handbar.
- For use with metal frame headboards only.
- Assisted weight capacity of 250 lbs.
- Designed to meet the requirements of HCPCS code: E0910.

2800A Chrome-plated 1 ea **2800GA** Gray 1 ea









