

Neb-u-Tyke™ Neb-a-Doodle

JOHN BUNN
A Graham-Field® Family Brand

This interactive, child-friendly design provides an advanced solution for effective aerosol treatment for children. Designed with an integrated drawing board, the Neb-a-Doodle will keep child relaxed and entertained during aerosol treatment and can be used over and over again!

The John Bunn Neb-u-Tyke Neb-a-Doodle is the latest addition to the Graham-Field line of quality nebulizer compressors that offer innovation, performance, and value for both patients and providers.



Features:

- Reliable, piston-driven compressor
- Interactive, drawing board for child to draw, doodle and write during aerosol treatment
- Convenient and compact design for easy transport and storage
- Comes complete with Nebulizer Accessory Kit
- Five year limited warranty
- Designed to meet the requirements of HCPCS Code: E0570

JB0112-070 1 ea

graham field



Neb-u-Tyke™ Neb-a-Doodle



- The John Bunn Neb-a-Doodle Nebulizer comes complete with:
- 1 Neb-u-Tyke Neb-a-Doodle Nebulizer Compressor
 - 1 User Manual
 - 1 Nebulizer Accessory Kit complete with:
 - 1 Reusable nebulizer (single patient use only)
 - 1 Adult mask
 - 1 Pediatric mask
 - 1 Mouthpiece
 - 1 PVC tubing
 - 1 "T" adapter
 - 3 Spare filters
 - 1 Reservoir tube

Technical Specifications

Power	115VAC, 60Hz
Power Consumption	≤ 60 watts
Operating Temperature Range	41°F TO 104°F
Operating Humidity Range	30 to 95% RH
Storage Temperature Range	-40°F TO 158°F
Storage Humidity Range	10 to 95% RH
Dimensions	10" x 6.5" x 4.2" (255 mm x 166mm x 108 mm)
Weight	4 lb
Max Pressure (Compressor)	32.5 psi
Max Flow (Compressor)	Approx 12 LPM
Particle Size (Nebulizer)	≤5 microns
Capacity (Nebulizer)	10 ml
Method Of Nebulization	Piston Driven
Sound Level	under 57 dBA

Information contained herein is subject to change.
The most current and complete product information can be found on our website.
Manufactured for GF Health Products, Inc. Made in China

Follow us on

graham field

www.grahamfield.com

© 2009, GF Health Products, Inc. All Rights Reserved.
John Bunn and Graham-Field are trademarks of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485 Certified Company.

The Graham-Field® Family of Brands
LUMEX Everest & Jennings Hausted Basic American JOHN BUNN

CE 0434

GF0900260RevB16