

The Everest & Jennings Bariatric Transport Chair is designed as an affordable bariatric patient transfer device and known for its extra width. With its 24" width and 400 lb weight capacity, this transport chair combines extra capacity in addition to compact transport need. Handle mounted attendant locks make it easy for a companion to engage.



Standard Features:

- Dual crossbrace with gusseted aluminum frame with epoxy painted finish in silver color
- Available in 24" width with 18" seat depth
- Black, padded nylon upholstery is flame retardant
- Fold-down back makes transport and storage easy
- Handle mounted attendant locks
- Large- flat free rear casters with 8" swivel front casters
- Removable flip-back desk length padded arms provide user comfort
- Detachable swingaway footrests for easy transfer
- Includes carry straps on seat for ease of lifting and folding
- Comes with positioning belt 42" length
- Available with anti-tippers on models after serial number IOS291
- Weighs 31 lb (w/o front riggings)
- Maximum Weight Capacity 400 lb – EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E1039



Warranty:

Three (3) year limited warranty

Bariatric Transport Chair		<i>Dimensions and Weights</i>									
		Seat			Back	Overall			Weight		
		Width	Depth	Height	Height	Depth	Height	Width	Folded Width	Without front rigging	Max Capacity
EJ777-3	Aluminum - Silver	24"	18"	20"	18"	40"	39 1/4"	30 1/4"	12"	31 lb	400 lb

EJ777-3 – Common Accessories & Parts		HCPCS Code
EJ76-BELT	Positioning belt (plastic buckle) 42" for Transport chairs	E0978
EJ777-3ANT	Anti-tippers, pair (for models after serial# IOS291)	E0971