

Everest & Jennings introduces the heavy duty articulating elevating legrests, industry compatible with aluminum footplate for durability. Normally used by tall individuals who require extension while fully elevated, this elevating legrest is sold as a separate accessory.



Standard Features:

- Aluminum footplate in black for bariatric or standard manual wheelchairs
- Extends in length from a minimum 16" in the downward position up to a maximum of 27" *when fully elevated*
- Swingaway detachable legrest can be adjusted for different needs or to avoid fatigue and enhance treatment
- Chrome-plated finish with black padded panel in vinyl for comfortable calf support
- Compatible with Advantage, Advantage Recliner, Traveler SE, Traveler LX, Traveler L3, Traveler L3 Plus, Traveler L4, Traveler HD & Traveler HTC
- Meets Medicare code K0053



Warranty:

One (1) year limited warranty

90763432 1 pr

LENGTH ADJUSTMENT RANGE*	
Position Downward	Fully Elevated Upward
Minimum 16"	Maximum 22"
Maximum 22"	Maximum 27"

*Measurements from seat upholstery edge to footplate