

# Everest & Jennings®

## WHEELCHAIRS

Setting standards in quality  
Enhanced transfer and comfort  
Manufactured for durability







# Everest & Jennings®

## Value:

Graham-Field® is known for a broad product offering at affordable pricing; without sacrifice to quality. Like you, we understand that the cost of a product is more than a price tag. The cost also includes the time and effort required to sell and service your customers. At Graham-Field we are proud to offer the value found in the Everest & Jennings™ line of wheelchairs.

## Quality:

E&J™ has maintained a commitment to provide our customers quality manual folding wheelchairs since 1933. In a time when many manufacturers are producing products at the lowest possible cost, often cutting corners on quality, E&J is standing by its dedication to quality through product design, specification control, testing and inspection. It is easy to see why the E&J Traveler® name is synonymous with quality.

## Service:

Our products are often only a day away with centralized distribution located in the metro areas of Atlanta, Georgia; St. Louis, Missouri; and Los Angeles, California. E&J also carries a full line of parts and upholstery to meet your service needs. In addition, our customer service department and technical service teams are always ready to help with your product needs.





# Manual Folding Wheelchairs

## Advantage®

The Advantage wheelchair provides comfort and durability at an affordable price. Designed as a low maintenance wheelchair.

### Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth
- Black padded nylon and vinyl upholstery available with chart pocket
- Carbon steel frame provides strength and durability
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Hammertone paint provides a long lasting and durable finish
- Comfortable padded armrests
- Maximum Weight Capacity – 250 lb – EVENLY DISTRIBUTED
- Meets Medicare code K0001
- One (1) year limited warranty



## Traveler® SE

The Traveler SE features a hammertone frame with saddle-joint tubular construction with dual axle for an adjustable 2" seat height and tool-free, industry compatible, front rigging.

### Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth and 20" by 18" seat depth
- Black vinyl upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Two-position axle and caster fork provide 2" seat height adjustment
- Maximum Weight Capacity 250 lb & 300 lb on 20"x18" models - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002 and K0006 on 20"x18" models
- One (1) year limited warranty



## Traveler® LX

The Traveler LX incorporates the legendary Everest & Jennings wheelchair chrome finish frame with saddle-joint tubular construction with dual axle for an adjustable 2" seat height and tool-free industry compatible front rigging. The premium vinyl upholstery is reinforced with an extra foam and heavy duty inner liner for enhanced comfort and durability.

### Standard Features

- Available in 16", 18" and 20" seat widths
- Blue premium vinyl upholstery with standard chart pocket on back
- Chrome finish provides a long lasting and easy to clean finish
- Tool-free front rigging with heel loops standard on swingaway footrests
- Gray non-marking tires
- Two-position axle and caster fork provide 2" seat height and adjustment
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002
- Two (2) year limited warranty





# Bariatric Folding Wheelchairs

## Traveler® HD

The Traveler HD, featuring a heavy-duty reinforced frame, provides a solution for users weighing up to 500 lb.

The Traveler HD comes with double crossbraces, heavy-duty 5/8" diameter axles, and a long lasting, embossed back and seat upholstery with a double inner liner for durability.

### Standard Features

- Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard)
- Embossed black vinyl upholstery with double inner lining for durability
- Hammertone paint provides a long lasting and durable finish
- Reinforced frame at all critical load-bearing points
- Chart pocket standard on upholstery
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007
- Five (5) year limited warranty



## Paramount XD

The Paramount XD is ideal for the bariatric wheelchair user. It features 26" and 30" seat widths as well as a weight capacity of 650 lb. The Paramount XD is constructed of 1" tubing to meet high weight capacity demands.

### Standard Features

- Available in 26" and 30" seat widths by 19 1/2" seat depth
- Swingaway footrest and elevating legrest options
- Maximum weight capacity of 650 lb - EVENLY DISTRIBUTED
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Reinforced welds and gussets on frame
- 8" x 2" Solid casters
- 24" Solid rear wheels
- Meets Medicare code K0007
- Five (5) year limited warranty



## Traveler® HTC

The Traveler HTC (Hospital Transport Chair) provides an alternative solution to the hospital/clinical transport need by incorporating special attendant and inventory control features with a heavy duty frame reinforced at all critical load-bearing points.

The Traveler HTC comes with an attendant push handle and rear mounted brakes that offer easier transport for the caregiver. In addition, the full width push handle acts as an anti-fold bar.

### Standard Features

- Available in 20", 22" and 24" seat widths, 18" seat depth (1" longer than standard)
- Fixed attendant push handle also functions as an anti-fold bar
- Includes a fixed anti-theft post ideal in hospital/clinical settings
- Rear mounted attendant wheel locks
- Fixed desk arms and front rigging lock nuts ensure inventory control
- Anti-tippers standard
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Black vinyl upholstery with chart pocket on back
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007
- Five (5) year limited warranty





# Specialty Folding Wheelchairs

---

## Advantage® Recliner

The Advantage Recliner comes with a higher back and removable full support headrest. The headrest includes an additional neck support that adjusts easily to fit every need. The durable hammertone frame is designed to support up to 300 lb and in the heavy duty version up to 450 lb users.

### Standard Features

- Available in 16", 18" and 20" seat widths, with 17" seat depth and in 22" seat width with 19" seat depth
- Features manual full reclining back range from 90-170 degrees
- Removable full support headrest with detachable padded desk or full arms
- Hammertone finish provides durability and style
- Embossed black leatherette type vinyl upholstery for easy cleaning
- Anti-tippers standard
- Comes with spreader metallic back bar which works as anti-fold device
- Industry compatible tool-free adjustable elevating front rigging
- Standard chart pocket on back
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED for 16", 18" and 20" models
- Maximum weight capacity 450 lb - EVENLY DISTRIBUTED for 22" models
- 16", 18", 20" models meet Medicare code K0001/ E1226 + E0966 + E0971 + K0195
- 22" models meet Medicare code K0007/E1226 + E0966 + E0971 + K0195 + E2201
- Five (5) year limited warranty



## Rehab Shower Commode

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or used as a portable, self-contained commode.

### Standard Features

- Available in 14 1/2" and 17" High back and Low back, caster pin locks and 5" rear locking caster models
- Completely removable full length arms
- Swingaway, detachable footrests
- Non-corrosive composite footplates
- Recessed seat and frame side contours
- Removable, hygienic plastic pail and lid standard
- Oval opening padded seat comes as standard
- Cushioned vinyl seat is sealed against moisture penetration
- Maximum weight capacity 250 lb – EVENLY DISTRIBUTED
- Meets Medicare code E0240
- Two (2) year limited warranty



**1.770.368.4700**

Information contained herein is subject to change.  
The most current and complete product information can be found on our website.  
**[www.grahamfield.com](http://www.grahamfield.com)**



# Lightweight Folding Wheelchairs

## Traveler® L3

The Traveler L3 lightweight frame weighs less than 36 lb making it the ideal wheelchair built for short and long term use. Designed with dual axles and multi-position caster forks, this wheelchair provides an adjustable seat height of 18" or 20". The wheelchair comes standard with non-marking polyurethane tires and industry compatible front rigging.



### Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth
- Black nylon upholstery with standard chart pocket on back
- Non-marking polyurethane rear tires
- Industry compatible tool-free front rigging with heel loops standard on swingaway footrests
- Two-position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003
- Five (5) year limited warranty

## Traveler® L3 Plus

The Traveler L3 Plus is an improved version of the previous Traveler L3 series. It is a lightweight Medicare approved K0003 wheelchair built for short or long term use. The wheelchair comes standard with non-marking polyurethane tires, industry compatible front rigging and adjustable anti-tippers along with a fully detachable 2" seat extension with hardware.



### Standard Features

- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Two position axle and caster fork provides 2" seat height adjustment
- Weighs less than 36 lb without front rigging
- Hammertone finish and carbon steel frame provide strength and durability
- Black nylon upholstery with standard chart pocket on back
- Velcro tab under seat to avoid gap between seat and back
- Anti-tippers come standard
- Tool-free adjustable front rigging with non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Fixed flip-back padded desk arms
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003
- Five (5) year limited warranty

## Traveler® L4

The Traveler L4 is a high-strength lightweight K0004 Medicare coded wheelchair. The Traveler L4 has dual axle and multi-position caster forks for an adjustable seat height of 18 or 20 inches. The wheelchair comes standard with non-marking polyurethane tires, anti-tippers and a chart pocket.

### Standard Features

- Available in 16", 18" and 20" seat widths, 16" seat depth. Also 18" and 20" seat widths with 18" seat depth
- Quick release axle is available for some models
- Black nylon upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Anti-tippers come standard
- Flip-back desk and full length arms to ease lateral transfer
- Height adjustable desk arms on some models
- Industry compatible tool-free front rigging with standard heel loops on swingaway footrests
- Two-position axle and caster fork provides 2" seat height adjustment
- Maximum weight capacity 250 lb EVENLY DISTRIBUTED (300 lb on models with 18" seat depth)
- Meets Medicare code K0004
- Five (5) year limited warranty





# Transport Chairs

## Aluminum Transport

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device and known for its durability. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

Rear mounted attendant locks make it easy for a companion to engage.

### Standard Features

- 17" and 19" seat widths
- Lightweight aluminum frame available in stylish color options of Aqua, Black, Red, Blue, Silver and Pink
- Fold down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Footrests are industry compatible and offer easy tool-free adjustment
- Full length padded arms provide user comfort
- Weighs 20 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code E10385
- One (1) year limited warranty



## Steel Transport

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device and known for its durability. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

Rear mounted attendant locks make it easy for a companion to engage.

### Standard Features

- 17" and 19" seat widths
- Steel frame is finished in a durable hammertone coating
- Fold down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Footrests are industry compatible and offer easy tool-free adjustment
- Full length padded arms provide user comfort
- Weighs 24 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Meets Medicare code E10385
- One (1) year limited warranty



## Bariatric Transport

The Everest & Jennings Bariatric Transport Chair is designed as an affordable bariatric patient transfer device. With its 24" width and 400 lb weight capacity, this chair satisfies additional compact transport needs with extra capacity. Handle mounted attendant locks make it easy for a companion to engage.

### Standard Features

- 24" seat width
- Lightweight aluminum frame with silver epoxy painted silver color finish
- Fold down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Standard Carry straps on seat allow for ease of lifting and folding
- Padded flip-back desk arms provide user comfort
- Weighs 31 lb (w/o front rigging)
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Meets Medicare code E1039
- Three (3) year limited warranty





# Accessories

## Anti-tippers



<b>90763565</b>	Anti-tipper, Advantage
<b>90763065</b>	Anti-tipper, Traveler SE & HD
<b>90763065</b>	Anti-tipper, Traveler L3
<b>90763230A</b>	Anti-tipper, Traveler L3 Plus

## IV/O2 Holders



<b>B10320F</b>	O2/IV Holder Combination (Advantage, SE, L3, L4, HD, LX, Paramount XD)
<b>B10325F</b>	O2 Holder (Advantage, SE, L3, L4, HD, LX, Paramount XD)

## Positioning Belts



<b>92133040</b>	Positioning Belt, Metal Buckle
<b>300-505</b>	Positioning Belt, Auto Style
<b>EJ76-BELT</b>	Positioning Belt, Plastic Buckle

## O2 Cylinder Holder Bag



<b>EJ02BAG07</b>	O2 Bag
------------------	--------

## Footrest



<b>90763030</b>	Footrest with Composite Footplate
<b>90763430</b>	Footrest with Aluminum Footplate

## High Pivot Legrest



<b>90763435</b>	High Pivot Elevating Legrest with Aluminum Footplate, Hammertone Finish
<b>90763135</b>	High Pivot Elevating Legrest with Aluminum Footplate, Chrome Finish

## Legrest



<b>90763031</b>	Elevating Legrest with Composite Footplate
<b>90763431</b>	Elevating Legrest with Aluminum Footplate

## Padded Seat



<b>90635070</b>	Oval Cutout
<b>90635175</b>	Open Front (C-shaped)

Item No.	Description
<b>Advantage</b>	
3H010220 / 3H010230	Advantage 16X16 Desk Arms, Footrest / Legrest
3H010120 / 3H010130	Advantage 18X16 Desk Arms, Footrest / Legrest
3H010320 / 3H010330	Advantage 20X16 Desk Arms, Footrest / Legrest
<b>Traveler SE</b>	
3E010220 / 3E010230	Traveler SE 16X16 Desk Arms, Footrest / Legrest
3E010120 / 3E010130	Traveler SE 18X16 Desk Arms, Footrest / Legrest
3E010320 / 3E010330	Traveler SE 20X16 Desk Arms, Footrest / Legrest
<b>Traveler LX</b>	
3E020120 / 3E020130	Traveler LX 18X16 Desk Arms, Footrest / Legrest
<b>Traveler L3</b>	
3F010220 / 3F010230	Traveler L3 16X16 Desk Arms, Footrest / Legrest
3F010120 / 3F010130	Traveler L3 18X16 Desk Arms, Footrest / Legrest
3F010320 / 3F010330	Traveler L3 20X16 Desk Arms, Footrest / Legrest
<b>Traveler L3 Plus</b>	
3F020220 / 3F020230	Traveler L3 PLUS 16X16 Desk Arms, Footrest / Legrest
3F020120 / 3F020130	Traveler L3 PLUS 18X16 Desk Arms, Footrest / Legrest
3F020320 / 3F020330	Traveler L3 PLUS 20X16 Desk Arms, Footrest / Legrest
<b>Traveler L4 (K0004)</b>	
3F020220 / 3F020230	Traveler L4 16X16 Desk Arms, Footrest / Legrest
3F020120 / 3F020130	Traveler L4 18X16 Desk Arms, Footrest / Legrest
3F020320 / 3F020330	Traveler L4 20X16 Desk Arms, Footrest / Legrest
<b>Traveler HD Heavy Duty</b>	
3G010320 / 3G010330	Traveler HD 20X18 Desk Arms, Footrest / Legrest
3G010420 / 3G010430	Traveler HD 22X18 Desk Arms, Footrest / Legrest
3G010520 / 3G010530	Traveler HD 24X18 Desk Arms, Footrest / Legrest
<b>Traveler Hospital Transport Chair (HTC)</b>	
3J010330	Traveler HTC 20X18 Desk Arms, Elevating Legrest
3J010430	Traveler HTC 22X18 Desk Arms, Elevating Legrest
3J010530	Traveler HTC 24X18 Desk Arms, Elevating Legrest
<b>Paramount XD Heavy Duty</b>	
5PX10620 / 5PX10630	Paramount XD 26W Desk Arms, Footrest / Legrest
5PX10820 / 5PX10830	Paramount XD 30W Desk Arms, Footrest / Legrest
<b>Aluminum Transport Chairs</b>	
EJ789-1 / EJ780-1	Transport Chair Green 17" / 19"
EJ781-1 / EJ782-1	Transport Chair Red 17" / 19"
EJ783-1 / EJ784-1	Transport Chair Silver 17" / 19"
EJ785-1 / EJ786-1	Transport Chair Black 17" / 19"
EJ787-1 / EJ788-1	Transport Chair Blue 17" / 19"
EJ790-1 / EJ791-1	Transport Chair Pink 17" / 19"
EJ793-1	Transport Chair Aqua 19"
<b>Steel Transport Chair</b>	
EJ796-1	Transport Chair Silver Vein Steel 19"
<b>Bariatric Transport</b>	
EJ777-3	Transport Bariatric Chair 24W
<b>Advantage Reclining Wheelchair</b>	
3K010230 / 3K010250	Advantage Recliner 16X17 Legrest, Desk Arms / Full Arms
3K010130 / 3K010150	Advantage Recliner 18X17 Legrest, Desk Arms / Full Arms
3K010330 / 3K010350	Advantage Recliner 20X17 Legrest, Desk Arms / Full Arms
3K020430 / 3K020450	Advantage Recliner 22x19 Legrest, Desk Arms / Full Arms
<b>Rehab Shower Commode</b>	
12022010 / 12022011	14 1/2" LoBack / 14 1/2" LoBack, Caster Pin Locks
12022013	14 1/2" LoBack with 5" Rear Locking Caster
12022110 / 12022111	17 1/2" HiBack / 17 1/2" HiBack, Caster Pin Locks
<b>*For a complete list of sizes and options available, please visit <a href="http://www.grahamfield.com">www.grahamfield.com</a></b>	