

The Traveler L3 Plus is an improved version of the previous Traveler L3 series. It is a lightweight Medicare approved K0003 wheelchair built for short or long term use. The wheelchair comes standard with non-marking polyurethane tires, industry compatible front rigging and adjustable anti-tippers along with a detachable 2" seat extension with hardware. Weighing less than 36 lb without front riggings, it features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers fixed flip-back desk length arm options for easy lateral movement.



2" seat extension with hardware comes standard

Standard Features

- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Two position axle and caster fork provides 2" seat height adjustment
- Weighs less than 36 lb without front rigging
- Hammertone finish and carbon steel frame provide strength and durability
- Black nylon upholstery with standard chart pocket on back
- Includes fastener tab under seat connecting both seat and back upholstery to eliminate gap
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Fixed flip-back padded desk arms
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003



Warranty

Five (5) year limited warranty













Chart pocket



Industry compatible front rigging



Anti-tippers come standard

Traveler® L3 Plus		Dimensions and Weights									
		Seat			Back	Overall				Weight	
		Width	Depth	Height	Height	Depth	Height	Width	Folded Width	Without front rigging	Max Capacity
											
3F012120	Flip-Back Desk Arm, Swingaway detachable footrest	18"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
3F012130	Flip-Back Desk Arm, Elevating legrest	18"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
3F012220	Flip-Back Desk Arm, Swingaway detachable footrest	16"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	24"	12 ¼"	34 lb	300 lb
3F012230	Flip-Back Desk Arm, Elevating legrest	16"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	24"	12 ¼"	34 lb	300 lb
3F012320	Flip-Back Desk Arm, Swingaway detachable footrest	20"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb
3F012330	Flip-Back Desk Arm, Elevating legrest	20"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb

Traveler® L3 Plus – Common Accessories & Parts		HCPCS Code
90763064A	Heel loop, pair	E0951
90763062A	Wheel lock lever extension, pair	E0961
90763SE16	2" seat extension kit for 16" W	E0981
90763SE18	2" seat extension kit for 18" W	E0981
90763SE20	2" seat extension kit for 20" W	E0981
300-505	Positioning belt 27-47" L (auto-style) push button	E0978
92133040	Positioning belt 36" L (airline style) metal buckle	E0978
90763230A	Anti-tippers adjustable, pair	E0971
907633CSE	Seat rail plug (dust cap) for seat extension, two pairs per bag	
B10320F	O2-IV holder	E2208
B10325	O2 cylinder holder	E2208
EJO2BAG07	O2 holder bag	E2208