

# Everest & Jennings<sup>®</sup> by graham-field

The surrounding positive effect in the home or facility can increase independence and allow patients to perform basic activities that we normally take for granted.

## From the Hospital to Home

Once a patient is discharged and sent home, there are needs that may require the use of a specific wheelchair. Everest & Jennings has you covered with various options to satisfy those needs for different users.



graham  
field<sup>™</sup>

Rehabilitation centers and the Veterans Administration (VA) have long worked to establish programs that specialize in the rehabilitation of SCI (Spinal Cord Injury) patients.

Learning to manage the everyday activities and care of patients with these injuries has become a growing priority and a way to improve their quality of life. This popular wheelchair has stood the test of time being an all-encompassing solution to many for over 30 years.



## Everest & Jennings® Rehab Shower Commode 1202 Series

- Designed to be used as a shower, commode and toilet chair all-in-one
- Made of corrosion-resistant stainless steel for durability and sturdiness
- It can be positioned over a standard toilet or be used as a portable, self-contained commode if the user prefers to self-propel
- Available in durable 24" rear wheels with handrims and solid rubber tires, with or without locking pin casters or with four 5" casters
- Low-back and Hi-back versions with completely removable full length arms for ease of transfer
- Available in Low-back 14¼", 16" and Hi-back 17¼", 19½"
- Available in 17" seat width and the NEW wider version 20" seat width
- Sealed vinyl cushioned oval seat and frame side contours allow easy access and comfort in toileting activities and bowel management procedures
- Removable, hygienic plastic pail and lid are included
- Non-corrosive composite footplates to resist water
- Maximum weight capacity **EVENLY DISTRIBUTED** is 250 lb on 14¼" and 17¼" back models and 300 lb on 16" and 19½" back models
- Meets Medicare Code E0240 specifications

**Everest & Jennings®**  
by graham-field

**graham**  
**field™**

**1.770.368.4700**

Information contained herein is subject to change.  
The most current and complete product information can be found on our website.

**[www.grahamfield.com](http://www.grahamfield.com)**

© 2018, GF Health Products, Inc. All Rights Reserved.  
Everest & Jennings is a trademark of GF Health Products, Inc.

Follow us on   

GF1800040RevA18