

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device. Its compact size and lighter weight allow a companion to lift and store the chair with ease.



Standard Features:

- Available in 19" seat width with black nylon upholstery
- Aluminum frame in stylish color options of Black and Blue Celeste
- Handle mounted attendant locks with cable
- 12" flat free rear wheels with 8" swivel front casters
- Fold-down back makes transport and storage easy
- Black padded nylon upholstery is flame retardant
- Fixed full length padded arms provide user comfort
- Tool-free adjustable detachable swingaway footrests with composite footplates
- Includes carry straps on seat for ease of lifting and folding and positioning belt 42" length
- Comes standard with 3" webbing black nylon leg strap with adjustable loops to slide over installed footrests
- Attendant wheel locks on both rear casters for better control and access by companion
- Weighs 21.38 lb (w/o front rigging)
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED



Warranty:

One (1) year limited warranty



Deluxe 12" Rear Wheel Aluminum Transport Chair		Dimensions and Weights									
		Seat			Back	Overall				Weight	
		Width	Depth	Height	Height	Depth	Height	Width	Folded Width	Without front rigging	Max Capacity
EJ870-1	Black	19"	16"	18"	18"	37"	38"	23"	10"	21.38 lb	300 lb
EJ871-1	Blue Celeste	19"	16"	18"	18"	37"	38"	23"	10"	21.38 lb	300 lb

Traveler® SE – Common Accessories & Parts		HCPCS Code
EJ76-BELT	Positioning belt (plastic buckle) 42" for Transport chairs	E0978