


The Lumex GF7200 5-Leg IV Stand sit on a 19" diameter, texturized, solid steel gray base that provides strength and durability. The IV stand travels on five 2" swivel steel housed casters for ease of maneuverability. It offers an aluminum post that has a height adjustment from 47" to 82" and is secured with a twist-lock adjustment collar. Available with two-hook or four-hook ram's horn options, this IV stand is designed to satisfy the needs for any home and facility.

**Features:**

- Solid steel, texturized 5 leg base provides strength, durability and superior support
- Height adjustable aluminum post from 47" – 82" that is secured with twist-lock adjustment collar
- Travels on five 2" swivel steel housed casters for ease of maneuverability
- Available with two-hook or four-hook ram's horn options
- Two (2) Year Limited Warranty



Locking Collar



Full Height



Low Height

Item	Description	UOM
GF72002	Two-Hook 5 Leg IV Stand	1 ea
GF72004	Four-Hook 5 Leg IV Stand	1 ea
GF7200P	IV Pole Only	1 ea
GF7200-LOCK	Replacement Locking Collet	4 ea/pk
GF7200-CST	Replacement Caster	1 pr



4-hook



2-hook