



## RAISED TOILET SEATS

# DELUXE TOILET SEAT RISERS

Essential item to assist with sitting and standing from a toilet seat.

### Features

- Adds 3.5" to toilet seat height, reducing bending effort
- Securely attaches to toilet bowl using bolts to prevent shifting
- Designed to work with your existing toilet seat and lid to better blend in with decor
- Available in round or elongated toilet bowl sizes

### With Arms (6482R-2, 6482ER-2)

- Push buttons allow for armrest to be added or removed easily

### Flip-up Seat (6481R-2, 6481ER-2)

- Designed to allow toilet seat riser to flip-up for easy cleaning



	overall Width		Seat Depth		Seat Height		Width Between Arms		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
6482R-2	22.5	57.15	17.25	43.82	3.5	8.89	19.5	49.53	300	136.08	2/cs
6482ER-2	23.5	59.69	19	48.26	3.5	8.89	20.5	52.07	300	136.08	2/cs
6481R-2	13.8	35.05	17.5	44.45	3.5	8.89	N/A	N/A	300	136.08	1 ea
6481ER-2	13.8	35.05	19	48.26	3.5	8.89	N/A	N/A	300	136.08	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2011, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.