



CLINICAL CARE RECLINERS

DELUXE CLINICAL CARE RECLINER

Three reclining positions for patient comfort. Plus a caregiver-activated Trendelenburg (or Shock) Position, is easily actuated from either side by a conveniently located foot pedal.

The wide, contoured seat and back, along with an integrated footrest, provide an enhanced level of patient comfort.

Wall saver armrest caps increase durability of the armrests and protect walls from marring.

Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster.

Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Side table surface is smooth and features tapered edges for easy cleaning.

Quick ship available.

Wide selection of fabrics and colors. See the [Swatch Card](#) for details.



Upholstered inside panels

Smooth tables with no-seam tapered edge

Options:

- Heat and Massage
- Heat Only feature
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow
- Drop-Arm (FR566DG)



Easily removed side panels for cleaning and maintenance



Trendelenburg Position can be actuated from either side of the recliner by a convenient foot pedal



Multiple Positions for Enhanced Patient Comfort

Model (xxxx represents fabric colors)	Wall Saver Armrest Edges	Fully Upholstered Arms	Heat & Massage	Heat Only	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Width between side tables in upright position		Depth fully reclined (Trendelenburg)		Weight Capacity EVENLY DISTRIBUTED		UOM
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
FR566G	●				31.25	79.4	37.5	95.3	49.5	125.7	20.5	52.07	19	48.3	22.25	56.5	8	20.3	21.75	55.2	51.5	130.8	75	190.5	400	181.44	1 ea
FR566GH	●	●			31.25	79.4	37.5	95.3	49.5	125.7	20.5	52.07	19	48.3	22.25	56.5	8	20.3	21.75	55.2	51.5	130.8	75	190.5	400	181.44	1 ea
FR566GH0	●		●		31.25	79.4	37.5	95.3	49.5	125.7	20.5	52.07	19	48.3	22.25	56.5	8	20.3	21.75	55.2	51.5	130.8	75	190.5	400	181.44	1 ea
FR566GxxxxPA	●			●	31.25	79.4	37.5	95.3	49.5	125.7	20.5	52.07	19	48.3	22.25	56.5	8	20.3	21.75	55.2	51.5	130.8	75	190.5	400	181.44	1 ea
FR566GHxxxxPA	●	●			31.25	79.4	37.5	95.3	49.5	125.7	20.5	52.07	19	48.3	22.25	56.5	8	20.3	21.75	55.2	51.5	130.8	75	190.5	400	181.44	1 ea
FR566GH0xxxxPA	●		●		31.25	79.4	37.5	95.3	49.5	125.7	20.5	52.07	19	48.3	22.25	56.5	8	20.3	21.75	55.2	51.5	130.8	75	190.5	400	181.44	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.