



CLINICAL CARE RECLINERS

POWERED BARIATRIC RECLINER

The Lumex Powered Bariatric Recliner is the perfect solution for larger patients who require a wider recliner with an increased weight capacity. This recliner, with a heavy-gauge steel frame, has a maximum weight capacity of 700 lb when evenly distributed! A four-function hand pendant activates the motor, allowing caregivers to safely recline and reposition patients.

Features

- Quick-release trigger mechanism on push handle allows recliner to swiftly move into Trendelenburg position
- Equipped with rechargeable battery pack for operation if the power goes out
- Four swivel casters facilitate transport and maneuverability in confined areas
- Durable, high impact side tables sit flush with chair arm to conveniently hold prep materials. Side table is smooth and features tapered edges for easy cleaning
- Legrest and backrest raise and lower independently, giving the caregiver complete control over positioning
- Comes complete with head pillow for extra support and patient comfort



Blue Ridge (427)

Doe Skin (851)

Jade (857)

Rosewood (863)



Multiple Positions for Enhanced Patient Comfort

| Model (xxxx represents fabric color) | overall Width | | overall Depth | | overall Height | | seat Width | | seat Depth | | seat Height | | Arm Height from seat | | Width between arms | | Width between side tables in upright position | | Depth fully reclined (Trendelenburg) | | Weight Capacity EVENLY DISTRIBUTED | | UOM |
|--|------------------|------|------------------|-------|-------------------|--------|---------------|-------|---------------|-------|----------------|-------|----------------------------|-------|-----------------------|-------|---|--------|--|-------|--|--------|------|
| | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | Kg | |
| FR588Wxxx | 38 | 96.5 | 59 | 149.9 | 54 | 137.16 | 26 | 66.04 | 24 | 60.96 | 27 | 68.58 | 8 | 20.32 | 24.4 | 61.98 | 48 | 121.92 | 77 | 195.6 | 700 | 317.51 | 1 ea |



—GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.