

FOLDING WALKERS


TWO BUTTON FOLDING WALKER WITH 5" WHEELS

Features

- Easy push dual-release folding mechanism enables user to fold walker sides independently
- Sturdy 1" diameter steel construction ensures maximum strength
- PVC contoured hand grip
- Comes with 5" wheels (Replacement leg set with wheels item #716280-LEGS)
- Rear glide tips allow for walker to move easily and smoothly over most surfaces
- Designed to meet the requirements of HCPCS code: E0143
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- One (1) year limited warranty

Replacement Parts

		UOM
716280-LEGS	Replacement Leg Set	1 set
716280-TIPS	Replacement Tips	1 pr



		Width inside handgrips		overall Width		Depth front to rear		minimum Height		maximum Height		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
716280A-4	Two Button Folding Walker	17.5	44.45	24.5	62.23	20.25	51.44	32.5	82.55	39.5	100.33	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-1900084Rev C22