

## MOBILITY ACCESSORIES

# AUTO-STOP™ WHEELS

### Features

- Designed as a wheel and braking system, Autostop wheels combine rolling and braking in one product
- For use on front legs of walkers, 5" rubber wheels are ideal for navigating difficult terrain such as grass, sidewalks and carpet. Wheels stop rolling when weight is applied to the front of the walker
- Recommended for use on the rear legs of the walker when fixed or swivel wheels are used on front walker legs
- Designed to meet the requirements of HCPCS code: E0155
- Attractive retail packaging
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty
- Kits of AutoStop Wheels and 3" or 5" fixed wheels also available

*These wheels are recommended for use on Lumex Walkers. (We cannot guarantee the performance of these accessories on other industry products.)*



	Wheel diameter		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	lb	Kg	
603450A	5	12.70	300	136.08	1 pr



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF2100093RevC22