

MOBILITY ACCESSORIES

TALL LEG EXTENSIONS

Features

- Increase the walker height by 4" (10.2 cm) for taller patients
- Adjustable in height for broad range of user heights
- Designed to meet the requirements of HCPCS code: E0158
- Designed to fit 1" standard walker tubing
- Limited lifetime warranty

These leg extensions are recommended for use on Lumex Walkers. (We cannot guarantee the performance of these accessories on other industry products.)



		UOM
603710A	Tall Leg Extensions	4 ea/set



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field