

SHOWER TRANSPORT CHAIRS

PVC RECLINING SHOWER CHAIR

Features

- Constructed of the highest quality health care grade PVC
- Threaded stem casters for additional safety
- Ergonomic curved rails provide support without sharp edges and reclining chairs have anti-slip handgrips for added safety
- Deluxe elongated extra heavy-duty, contoured open front soft seat
- Padded lumbar support and head section.
- Dual swing away arms for patient transfer can be secured with red safety straps
- Elevated Leg Rest (89440 and 89550 only)
- Mesh available in your choice of 3 colors -
- Royal Blue mesh will be provided for online orders. If you desire another color, please order via Customer Service at 770-368-4700
- One (1) year limited warranty on PVC frame; six (6) month limited warranty on all other materials



Standard Mesh and Solid Vinyl Colors



Royal Blue



Forest Green



Mauve

	internal Width		external Width		Height		Depth		Seat Height		Back Width exterior		Length Upright (leg down)	Length Reclined (leg up)	Reclined Back		Weight Capacity EVENLY DISTRIBUTED		UOM		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg			
89330	20	50.8	24	60.96	48	121.92	18	45.72	21	53.34	16	40.64	43	109.22	58	147.32	29	73.66	325	147.42	1 ea
89440	26	66	30	76.2	50	127	22	55.88	22	55.88	23.25	59.06	41	104.14	62	157.48	27	68.58	600	272.16	1 ea
89550	30	76.2	34	86.36	50	127	22	55.88	22	55.88	27.25	69.22	41	104.14	62	157.48	27	68.58	700	317.51	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0700129RevC22