



## PURETILT

The PureTilt® patented Tilt-in-Space Wheelchair, designed by an Occupational Therapist to provide safer and independent transitions to facilitate optimal comfort for a patient.

A manual tilt wheelchair that has a user-adjusting tilt control so that the wheelchair user can readjust their own posture throughout the day. This is vital for comfort and pressure relief.

### Features:

- Height-adjustable, ergonomic, contoured back with adjustable contoured headrest
- User / caregiver-operated seat-tilt adjustment controlled by a built-in lever that uses a gas spring piston
- Padded, hinged, non-detachable arms swing up and back to facilitate transfer
- Back folds forward to facilitate storage and transport
- 22-inch rear wheels with aluminum hand rims for easy self-propulsion
- Two front casters swivel for easy maneuvering
- Rear anti-tippers standard
- Angle-adjustable leg rest assembly with:
  - User-operated mechanism with telescoping adjustable length
  - Footplate folds flat and locks in place
- HCPCS code E1161
- 300 lb (136.08 kg) Weight Capacity EVENLY DISTRIBUTED



### Accessories: sold separately

- Item No. PT3000-SP — Right & Left Side Pads
- Item No. PT3000-FP — Larger Leg Pad & Foot Pad

side pad set and foot pad accessories shown left - each sold separately

### Specifications:



folds for storage and vehicle transport

Item	Overall Chair Width		Overall Chair Length				Storage Dimensions with back folded				Seat Width		Seat Depth		Seat Height front of seat to floor				Main, Seat, and Back Frame	Wheels, Legrest, Footplate	Other Durable Components <sup>6</sup>
	inch	cm	upright		reclined		Length		Height with seat		inch	cm	inch	cm	bottom pin		top pin				
PT3000-18	28.50	72.39	39.80	101.09	50.00	127.00	39.50	100.33	36.00	91.44	19.80	50.29	21.00	53.34	19.75	50.17	21.75	55.25	60	12	6

  

Item	Arm Height		Foot Support Length distance from floor				Front Caster Diameter		Rear Wheel Diameter		Seat Tilt Leg to Seat <sup>5</sup> adjustment <sup>2</sup>		Shipping <sup>3</sup> Weight		Product Weight with cushions		Weight Capacity <small>EVENLY DISTRIBUTED</small>		warranty (months)
	inch	cm	retracted		extended		inch	cm	inch	cm	degree of movement		lb	kg	lb	kg	lb	kg	
PT3000-18	5.80	14.73	13.00	33.02	17.50	44.45	6.00	15.20	22.00	55.90	5 to 40	-97 to 15	100	45.36	95	43.09	300	136.08	

<sup>1</sup> ± 0.5 inches (± 1.27 cm) <sup>2</sup> ± 2° <sup>3</sup> estimated <sup>4</sup> not otherwise listed <sup>5</sup> from horizontal



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved. Everest & Jennings is a registered trademark of GF Health Products, Inc.