



FOAM POSITIONING PRODUCTS



BED WEDGES

The Lumex[®] Bed Wedges can be used to elevate legs, head or shoulder to help relieve pressure during patient recovery.

Features

- 100% polyurethane foam
- Removable and washable cover
- Latex Free



Specifications

Item	Description	Overall Width		Overall Length		Overall Height		Product ¹ Weight		UOM	Color	Foldable	Made in USA	Product Link	Warranty in months
		inch	cm	inch	cm	inch	cm	lb	kg						
LE247	Folding Bed Wedge	25	63.50	24	60.96	7	17.78	1.13	0.513	1/ea	Blue	■	—	Product Link	3
7-2394R	7 inch Bed Wedge	24	60.96	24	60.96	7	17.78	1.13	0.513	1/ea	White	—	■	Product Link	3
7-2393R	10 inch Bed Wedge	24	60.96	24	60.96	10	25.40	2.70	1.225	1/ea	White	—	■	Product Link	3
7-2395R	12 inch Bed Wedge	24	60.96	24	60.96	12	30.48	2.12	0.962	1/ea	White	—	■	Product Link	3

■ Standard ■ Optional ■ Available ■ Not Available ■ Accessory — Does Not Apply [Active Hyperlink](#)

¹ approximate product weight — as seat width increases, so does product weight, however slightly



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved. Lumex is a registered trademark of GF Health Products, Inc.