

LIFT SLINGS

HOYER COMPATIBLE PADDED SLINGS

Features

- Deluxe closed cell padding with four point loops for comfort and safety
- Available with or without head support
- Durable polyester material is manufactured to the same weight and size specifications as the name brand
- Value priced without compromising on quality
- Designed to meet the requirements of HCPCS code: E0621
- Six (6) month limited warranty



DSHC70012



DSHC70002

	Head Support	Sling Size	Weight Capacity EVENLY DISTRIBUTED		optimal User Weight		UOM
			lb	Kg	lb	Kg	
DSHC70000	—	XL	600	272.2	270 - 600	122.5 - 272.2	1 ea
DSHC70001	—	Large	500	226.8	198 - 350	89.8 - 158.8	1 ea
DSHC70002	—	Medium	500	226.8	99 - 210	44.9 - 95.3	1 ea
DSHC70010	■	XL	600	272.2	270 - 600	122.5 - 272.2	1 ea
DSHC70011	■	Large	500	226.8	198 - 350	89.8 - 158.8	1 ea
DSHC70012	■	Medium	500	226.8	99 - 210	44.9 - 95.3	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0900076RevB23