

LAMPS

E-SERIES HAND HELD UV-A LAMPS

Features

- Featuring a silver-anodized aluminum housing that makes the lamp both attractive and durable
- Full range of models offers a selection of UV intensity, size and wattage
- Perfect for applications where a compact lightweight lamp is desired



2212



2211



2213



| Tube | | Typical Peak UV Intensity at 6" (15.2 cm) | overall Length | | overall Width | | overall Height | | UOM |
|------|-------------------------------------|---|-------------------|------|------------------|-----|-------------------|-----|------|
| | | μW/cm ² | inch | mm | inch | mm | inch | mm | |
| 2211 | (1) 4-Watt integrally filtered tube | 800 | 7 | 17.8 | 3.25 | 8.3 | 2.5 | 6.4 | 1 ea |
| 2212 | (1) 6-Watt integrally filtered tube | 1,280 | 10.25 | 26. | 3.25 | 8.3 | 2.5 | 6.4 | 1 ea |
| 2213 | (2) 4-Watt integrally filtered tube | 620 | 7 | 17.8 | 3.25 | 8.3 | 2.5 | 6.4 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a registered trademark of GF Health Products, Inc.