

Cardiology Dual-Frequency Stethoscope

- Stainless steel single-sided chestpiece with dual-frequency diaphragm
- Specially designed membrane lock, within the chestpiece, alters the diaphragm to achieve dual-frequency
- Hold the chestpiece to achieve low frequency response (similar to dual-sided bell chestpiece)
- Apply firm pressure to achieve high-frequency diaphragm response
- Comfortable, non-chill ring
- Lightweight, stainless steel binaurals with internal spring
- Extra-thick cardiology 19" latex-free Y-tubing
- Deluxe eartips
- Lifetime Warranty

435 Black 1 ea



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0700061RevA07