

## Alternating Pressure/Low Air Loss Overlay System

## 5" Alternating Pressure/Low Air Loss Overlay System

- Slim, lightweight system offers a cost-effective alternative to alternating pressure mattress replacement systems
- 5" cells are constructed with soft yet durable Nylon/Polyurethane
- Nine laser cut low air loss cells in the sacral region are designed to keep patient cool and dry
- Pump adjustment from soft to firm increases patient comfort
- Includes low pressure indicator
- 265 lb maximum weight capacity
- Designed to meet the requirements of HCPCS code: E0372
- 1 year limited warranty
- Latex Free

14025 Alternating Pressure/Low Air Low Mattress Overlay and Pump System









