

# Everest & Jennings®



## Built to Last



GRAHAM-FIELD

Everest & Jennings strives to create wheelchairs that are comfortable, affordable and built to last. Offering various seat sizes, a reinforced frame and a 500 lb **EVENLY DISTRIBUTED** weight capacity has proven the Traveler HD to be a reliable option for your hospital needs.

Each wheelchair is carefully engineered and goes through rigorous testing to ensure stability and addresses the needs of both the patient and health-care professional.

**Embossed leatherette upholstery with standard chart pocket** include double inner lining to minimize sagging and maximize wear time.



**Silvervein, hammertone finish** is durable and more resistant to scratches upon impact.

**Padded and upholstered armrests** offer comfort and are available in both desk and full arm styles.

**Polyurethane tires** and composite handrims are long lasting and support heavy weight users.

**Dual crossbraces** with frame reinforcements at critical weight bearing points offer stability.

**Aluminum footplates** on footrest and elevating legrests are resilient with standard heel loops on footrest models.



**Anti-Fold Theft Bar**  
92523850

**O2 Holder**  
B10325



**Anti-Removal Fasteners**  
9B761770

**O2/IV Combo**  
B10320F



**Ask about our Custom Color Upholstery Program**

<b>Traveler® HD</b>		<i>Dimensions and Weights</i>								
		Seat			Back	Overall				Wt.*
<b>3G010320</b>	Detachable Desk Arm, Swingaway Footrest	20"	18"	18" or 20"	16 1/2"	34"	36"	28 1/2"	12 1/2"	57 lbs
<b>3G010330</b>	Detachable Desk Arm, Elevating Legrest	20"	18"	18" or 20"	16 1/2"	34"	36"	28 1/2"	12 1/2"	57 lbs
<b>3G010340</b>	Detachable Full Arm, Swingaway Footrest	20"	18"	18" or 20"	16 1/2"	34"	36"	28 1/2"	12 1/2"	57 lbs
<b>3G010350</b>	Detachable Full Arm, Elevating Legrest	20"	18"	18" or 20"	16 1/2"	34"	36"	28 1/2"	12 1/2"	57 lbs
<b>3G010420</b>	Detachable Desk Arm, Swingaway Footrest	22"	18"	18" or 20"	16 1/2"	34"	36"	30 1/2"	12 1/2"	57 lbs
<b>3G010430</b>	Detachable Desk Arm, Elevating Legrest	22"	18"	18" or 20"	16 1/2"	34"	36"	30 1/2"	12 1/2"	57 lbs
<b>3G010440</b>	Detachable Full Arm, Swingaway Footrest	22"	18"	18" or 20"	16 1/2"	34"	36"	30 1/2"	12 1/2"	57 lbs
<b>3G010450</b>	Detachable Full Arm, Elevating Legrest	22"	18"	18" or 20"	16 1/2"	34"	36"	30 1/2"	12 1/2"	57 lbs
<b>3G010520</b>	Detachable Desk Arm, Swingaway Footrest	24"	18"	18" or 20"	16 1/2"	34"	36"	32 1/2"	12 1/2"	57 lbs
<b>3G010530</b>	Detachable Desk Arm, Elevating Legrest	24"	18"	18" or 20"	16 1/2"	34"	36"	32 1/2"	12 1/2"	57 lbs
<b>3G010540</b>	Detachable Full Arm, Swingaway Footrest	24"	18"	18" or 20"	16 1/2"	34"	36"	32 1/2"	12 1/2"	57 lbs
<b>3G010550</b>	Detachable Full Arm, Elevating Legrest	24"	18"	18" or 20"	16 1/2"	34"	36"	32 1/2"	12 1/2"	57 lbs
<b>3G020500</b>	Fixed Full Arm, Fixed Footrest	24"	18"	18" or 20"	16 1/2"	34"	36"	32 1/2"	12 1/2"	57 lbs

Visit our website to see the full E&J Wheelchair line [www.grahamfield.com](http://www.grahamfield.com) or call 770.368.4700