




Traveler[®] HD

PARTS CATALOG



CONTENTS

SIDEFRAME AND CROSSBRACE3

ARM, DETACHABLE4

UPHOLSTERY4

REAR WHEEL5

CASTER.....5

FRONT RIGGING, DETACHABLE6

ACCESSORIES7

TO ORDER

Tel: 770-368-4700 • Fax: 770-368-2386


Returns: Parts may not be returned without prior authorization from GF Health Products, Inc. Custom parts cannot be returned for credit. Please refer to GF Health Products, Inc. Terms and Conditions for further information.

GF Health Products, Inc. is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this catalog are based on the latest product information available at the time of printing. The most current product information, including the most current revision of this catalog, can be found online at www.grahamfield.com.

Everest & Jennings, E&J, Traveler, and Graham-Field are registered trademarks of GF Health Products, Inc.

© 2009, GF Health Products, Inc., All Rights Reserved.

SIDEFRAME AND CROSSBRACE



1 Sideframe, all widths

Right 90763401
Left 90763402

2 Crossbrace, pair

Includes slides, center bolt and nut

20" W 907634C1
22" W 907634C2
24" W 907634C3

3 Bottom inner rail kit

Includes bottom inner rail (pair),
screw and nut that hold bottom
inner rail (4 each), and sleeve (pair) ... 90763503


4 Grip (10 each) 90763005

5 Tipping sleeve (10 each)..... 90763006

6 Seat rail cradle, screw, washer and
nut (pair)..... 90763403


7 Seat slide and seat rail sleeve
(2 each)..... 90763404

ARM, DETACHABLE



- 1 Arm, detachable fixed-height, pair
Includes right and left arm with side panel and spring lock (does not include arm pad)
Desk length..... 90763020
Full length..... 90763021
- 2 Arm pad, pair
Includes two black arm pads and mounting screws
Desk length..... 90763022
Full length..... 90763023
- 3 Spring lock (10 each) 90763024

UPHOLSTERY




- Back upholstery, black leatherette,
with built-in chart holder
Includes mounting strips and screws

20" W..... 907634A1
22" W..... 907634A2
24" W..... 907634A3

- Seat upholstery, black leatherette
Includes mounting strips and screws


20" W x 18" D..... 907634B1
22" W x 18" D..... 907634B2
24" W x 18" D..... 907634B3

REAR WHEEL




- 1 Rear wheel with handrim, 24" x 1",
polyurethane tire
Includes hub cap and bearing;
does not include axle or nut..... 90763450
- 2 Hub cap (10 each)..... 90763451
- 3 Wheel bearing (10 each) 90763452
- 4 Rear axle and nut (10 each)..... 90763453
- 5 Spacer (10 each)..... 90763454

CASTER




- Complete caster (fork and wheel),
8" x 2" solid urethane tire 90763440
- 1 Caster wheel, axle and nut 90763441
- 2 Caster fork, stem nut and dust cap 90763042
- 3 Caster stem bearing (10 each) 90763043
- 4 Caster dust cap (10 each)..... 90763044
- 5 Caster axle and nut (10 each)..... 90763045

FRONT RIGGING, DETACHABLE




Footrest, pair, all widths
Includes right and left hanger with release
and aluminum footplate with extension
and heel loop 90763430

Aluminum footplate, pair
(Fits all front rigging) 90763433



Low-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
aluminum footplate with extension, and
legrest panel (does not include
heel loop) 90763431

Elevating legrest panel kit
(Fits all elevating legrests)
Includes 2 black padded panels
and hardware 90763032



High-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
aluminum footplate with extension, and
legrest panel (does not include
heel loop) 90763435

Traveler HD accessories (not shown)

1	Heel loop, pair.....	90763064
2	Wheel lock, pair	
	Includes rubber tips	
	Push to lock.....	90763060
	Pull to lock.....	90763061
3	Rubber tip, wheel lock (10 each).....	90763063
4	Wheel lock extension, pair	90763062
5	Positioning belt (metal buckle)	90133125
6	Anti-Fold bar	92523850
7	O2-IV holder	B10320F
8	Anti-Tippers (pair)	
	20" seat height	90763065
	18" seat height	90763066

Everest & Jennings®

**GF Health Products, Inc.
2935 Northeast Parkway
Atlanta, Georgia 30360 USA**

tel: 770-368-4700

fax: 770-368-2386



GRAHAM-FIELD

www.grahamfield.com

**Traveler HD Parts Catalog
GF0900210RevB14**