

The Lumex Set n' Go Wide Height Adjustable Rollator offers the convenience and flexibility of a two-in-one product all rolled into one with an added offering for larger individuals. With its enhanced, wider spacing of 19" in between the handgrips and larger 15.75" seat, the Set n' Go Wide is ideal for those users needing a little more room for comfortable movement and seating. Designed with multiple height settings on the frame, the Set n' Go Wide allows the seat height to be adjusted from a hemi 19.5" to a standard 23.5" height in 1" increments. This accommodates for a broad range of user heights which in turn offers more comfort for the user when seated.

**2 Rollators in 1!**



**19" WIDE**

## Features:

- Easy height adjustment on the frame allows seat height to be adjusted from 19.5" to 23.5" in 1" increments
- Comfortable, 15.75" wide, seamless padded seat with zippered pouch under seat for security of personal items
- Also includes black pouch for storage of larger items
- Ergonomic handgrips offer a comfortable grip
- Adjustable handle height with knobs to secure in place
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Removable padded backbar with easy release pull knobs
- 8" wheel for both indoor and outdoor use
- 350 lb maximum weight capacity - EVENLY DISTRIBUTED
- Velcro strap attached secures rollator frame when folded for storage and transport
- Designed to meet the requirements of HCPCS: E0143 and E0156
- Five (5) year limited warranty



- RJ4718B** Blue
- RJ4718R** Burgundy
- RJ4718S** Silver



Folding Strap

**Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.**

## Specifications

### Width

Inside Handgrips: ..... 19.4" (492 mm)  
Overall: ..... 26.9" (683 mm)

### Depth

At Minimum Seat Height: ..... 27.5" (699 mm)  
At Maximum Seat Height: ..... 31.5" (800 mm)

### Height

At Minimum Handle Height: ..... 32" (813 mm)  
At Maximum Handle Height: ..... 38" (965 mm)

### Seat

Height: ..... 19.5" – 23.5" (495 mm - 597 mm)  
Width: ..... 15.75" (400 mm)  
Depth: ..... 12" (305 mm)  
Casters Diameter: ..... 8" (203 mm)  
Maximum User Weight: ..... 350 lb (159 kg)

EVENLY DISTRIBUTED

Product Weight: ..... 16.5 lb (7.5 kg)