

Fast Alert Patient Alarm with Magnetic Pull Cord

- Economical alarm that notifies staff or caregiver that a user has moved
- Alligator clip attaches to clothing so when user gets up, the cord pulls from the alarm and the audio alarm is activated
- Activation cord is adjustable from 17" to 31"
- Alarm easily secures with self contained clip
- Convenient dual volume settings: 70dB and 90dB
- Requires one 9V battery (not included)
- Low battery indicator
- 1 Year Limited Warranty

GF13700 1 ea

