

The Grafco 6-Place Fixed-Angle Laboratory Centrifuge is an affordable fixed-angle centrifuge designed for blood separation in a small lab or physician office. The compact design uses less than 1 square foot of counter space while delivering high-quality fixed angle separations. Heavy gauge, durable steel construction and reinforced side walls offer exceptional strength and safety protection for a durable and dependable unit that will last for years. Its maintenance-free brushless motor and custom tube holders provide reliable and consistent operation without requiring additional accessories. The centrifuge helps maintain a quiet laboratory setting, performing at just 55 dB.



- **Easy to Clean:** Removable rotor and accessories can be easily cleaned and disinfected with isopropyl alcohol, soap and water, or a 10% (5500 PPM) bleach solution
- **Cool Running:** Motor housing and rotation chamber allows for cool operation under the heaviest workloads
- **Standards & Certifications:** Conforms to the latest clinical and safety standards, including UL and CAN/CSA**
- **Clear Shatter-proof Lid:** Offers safe sample observation and optical speed calibration
- **Lid Safety Latch:** Unit will not operate unless lid is closed and latched.
In the event that the lid is opened during a run, the power to the motor is disconnected
- **Brushless AC motor:** Provides years of operation with no routine maintenance
- **Quiet operation:** Double-cased motor and stability mount design provide stability against imbalance
- **Mechanical Timer:** 1 to 30 minutes allows for wide range of applications
- **Capacity:** 90 ml (6 x 15ml test tubes up to 17mm x 125mm)
- **Overall Dimensions (H x W x D):** 8.75" x 11.75" x 14" (222.3 mm x 298.5 mm x 355.6 mm)
- **Power Requirement:** 115 V AC, 60 Hz
- **Weight:** 11 lb (4.99 kg)
- **Nominal Acceleration Time:** 45 seconds
- **Max Acceleration Time:** 30 seconds
- **Max Deceleration Time:** 40 seconds
- Rated for a rotation frequency of 3,300 RPM
- One (1) Year Limited Warranty



Includes:

One (1) six-place fixed-angle rotor centrifuge
Six (6) 75 mm tube holders (Item GF614B-75)
Six (6) 100 mm tube holders (Item GF614B-100)
Six (6) 125 mm tube holders (Item GF614B-125)
User Manual

**This product is not made with natural rubber latex*

Item #	Description	Color	UOM
GF614B	6-Place Fixed Angle Centrifuge		1 ea
Optional Accessories:			
GF614B-75	75 mm tube holders	Green	6/bx
GF614B-100	100 mm tube holders	Red	6/bx
GF614B-125	125 mm tube holders	Black	6/bx

Tube Holders	Nominal Speeds (RPM)	Nominal RCF (x g)
with 125 mm holders	3150 (+/- 100)	1200 (+/- 80)
with 100 mm holders	3250 (+/- 100)	1080 (+/- 80)
with 75 mm holders	3250 (+/- 100)	950 (+/- 80)

** Complies with all requirements of UL standard 61010A-1, 61010A-2-20; CAN/CSA C22.2 No's 1010.1; 1010.2.20.