

Automatic Gas Lift Overbed Table Base

The Acute Care overbed table features a state-of-the-art gas cylinder lift mechanism that provides quiet, smooth, and reliable movement and can be mounted on either the right or left side of the base. The unique lift mechanism provides for infinite height adjustment between 27" and 42" and allows the table to be adjusted upward without requiring the handle to be activated. This convenience provides extra security if the bed rises with the table over it.

Capable of handling 100 lb of evenly distributed weight. The heavy duty column and base come in a soft, neutral, Opal powder coat finish for aesthetic appeal that coordinates with almost any decor. The base is available in either an "H", "U", or "Fork" design. The oversized 2" dual wheel casters provide for easy maneuverability and movement over thresholds.

Item: **OB4828NT**



Double Mirror Vanity

Made of durable polymer construction, this dual compartment vanity can be accessed from either side of the table top. This model features two large, hinged, acrylic mirrors attached with a living hinge, and rounded corners for easy cleaning. Field installed on tabletop.



Item: **OB4828VAN**

Tabletops

Standard Overbed tabletop 15" x 30" A28
Oversized Overbed tabletop 18" x 35" A28A

Standard Decorative Laminates

