



**This product is not made with natural rubber latex*

- Medium acuity mattress for patients with low grade ulcers
- Convoluted 5-zone "Wave" foam top surface for increased airflow and comfort
- Tri-Laminate construction with soft heel-zone
- Firm bolster with optional 3" raised bolster
- Meets 16 CFR 1633
- Designed to meet the requirements of HCPCS Code E0184
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED
- Four (4) year limited warranty on mattress,
One (1) year limited warranty on cover

Dimensions	Flush Bolster	Raised Bolster
76" X 35"	C400B-3576	C400RB-3576
80" X 35"	C400B-3580	C400RB-3580
84" X 35"	C400B-3584	C400RB-3584
76" X 42"	C400B-4276	C400RB-4276
80" X 42"	C400B-4280	C400RB-4280
84" X 42"	C400B-4284	C400RB-4284