

The LUMEX Select Convertible Mattress provides outstanding pressure displacement and redistribution to assist in the prevention and treatment of stage I through IV pressure ulcers.

Features:

- Closed System (Static Air) can quickly convert to an Alternating Pressure Mattress
- Eight non-circular cylinders prevent "rolling" or "gapping"* and are fully enclosed in a fiber mesh sleeve and jacket
- Air control unit is simple, quiet and set at the optimum 10 minute cycle to eliminate guesswork for the caregiver
- 2.5" Hybrid Visco-elastic top layer does not raise body heat as with standard visco-elastic foam
- Stretch cover is fluid resistant, air-permeable (breathable) and reduces friction and patient shear

Use of the Lumex Select Convertible will eliminate the need for:

- costly rental and storage of mattresses that are only used on an as needed basis
- moving and transporting mattresses from room to room
- overlays that are low in profile and do not provide adequate support or therapeutic benefit



Specifications

Standard Size: 35" x 80" x 7"
 Custom sizes available including bariatric sizes (42" - 48")
 Maximum Weight Capacity 500 lb (750 lb on bariatric models) - EVENLY DISTRIBUTED
 Consult frame manufacturer regarding capacity
 Compliant with 1633, Cal TB 129, ASTM E 1590 and BFDIX-11 (Boston) fire safety regulations
 Available with 2" side-bolsters
 Five (5) year warranty on mattress

**This product is not made with natural rubber latex*

LCM3580	35" x 80", Convertible Non-Powered Mattress	1 ea
LCM3580RB	35" x 80", Convertible Non Powered Self Adjusting Convertible Mattress with Side Bolsters	1 ea
LCM4280	42" x 80", Convertible Non-Powered Mattress	1 ea
LCM4280RB	42" x 80", Convertible Non Powered Self Adjusting Convertible Mattress with Side Bolsters	1 ea
LCM4280B	42" x 80", Convertible Non-Powered Mattress-Bariatric	1 ea
LCM4880	48" x 80", Convertible Non Powered Self Adjusting Convertible Mattress	1 ea
LCM4880RB	48" x 80", Convertible Non Powered Self Adjusting Convertible Mattress with Side Bolsters	1 ea
LCM4880B	48" x 80", Convertible Non Powered Self Adjusting Convertible Mattress- Bariatric	1 ea
LCMP100	Pump	1 ea