

# LUMEX®

A Graham-Field® Family Brand

## Walkabout Essentials Four Wheel Rollator

The Lumex Walkabout Essentials Four Wheel rollator offers a comfortable and cost effective solution to mobility for someone looking for an easy to use walking aid. The frame is made of a highly durable steel with 6" easy to maneuver casters. The Walkabout quickly and easily folds into a compact unit for storage and transport.



### Features:

- Ergonomic hand grips make them easy to grip
- Adjustable handle height to accommodate broad range of user heights
- Durable steel frame
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable, hinged, padded backrest can be folded up and down
- 6" casters for ease of maneuverability
- Pouch included for storage of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Three (3) year limited warranty



**RJ5000B** Blue 1 ea  
**RJ5000R** Red 1 ea



Color options

**Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.**

### Specifications:

Width Inside Handgrips .....	15.5"
Overall Width .....	24"
Overall Depth .....	30.7"
Seat Width .....	13"
Seat Depth .....	12"
Seat Height .....	19.75"
Caster Diameter .....	6"
Minimum Handle Height .....	30"
Maximum Handle Height .....	36"
Product Weight .....	17.4 lb (7.9 kg)
Weight Capacity .....	300 lb

EVENLY DISTRIBUTED

**graham**  
**field**

[www.grahamfield.com](http://www.grahamfield.com)

© 2017 GF Health Products, Inc.  
Atlanta, GA All rights reserved.  
GF1700089RevA17