

Everest & Jennings introduces the Traveler® L3 Plus as the lightweight wheelchair addition to the Traveler® line, a brand respected for performance and durability. The Traveler® L3 Plus lightweight frame weighs less than 36 pounds making it the ideal K0003 Medicare coded wheelchair built for short and long term use as well as rental needs.

This lightweight wheelchair comes standard with non-marking polyurethane tires, industry compatible front rigging and adjustable anti-tippers along with a detachable 2" seat extension with hardware. It features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers flip-back desk, adjustable height and full arm options for easy lateral movement. The Traveler® L3 Plus is manufactured with the same commitment to quality Everest & Jennings is known for.



2" seat extension with hardware comes standard



Anti-tippers come standard

### Standard Features

- Weighs less than 36 lb without front rigging
- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Two position axle and caster fork provides 2" seat height adjustment
- Hammetone finish and carbon steel frame provide strength and durability
- Black nylon upholstery with standard chart pocket on back
- Includes fastener tab under seat connecting both seat and back upholstery to eliminate gap
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Flip-back padded arms available in fixed desk length, adjustable height desk and detachable flip-back full length
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003



### Warranty

Five (5) year limited warranty













Chart pocket



Full length arms



Adjustable height arms

<b>Traveler® L3 Plus</b>		<i>Dimensions and Weights</i>									
		<i>Seat</i>			<i>Back</i>	<i>Overall</i>				<i>Weight</i>	
		<i>Width</i>	<i>Depth</i>	<i>Height</i>	<i>Height</i>	<i>Depth</i>	<i>Height</i>	<i>Width</i>	<i>Folded Width</i>	<i>Without front rigging</i>	<i>Max Capacity</i>
											
<b>3F012120</b>	Flip-Back Desk Arm, Swingaway detachable footrest	18"	16"*	18" or 20"	17 ¼"	29"	35½"	26"	12 ¼"	34 lb	300 lb
<b>3F012130</b>	Flip-Back Desk Arm, Elevating legrest	18"	16"*	18" or 20"	17 ¼"	29"	35½"	26"	12 ¼"	34 lb	300 lb
<b>3F012140</b>	Flip-back Full Arm, Swingaway Footrest	18"	16"*	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
<b>3F012220</b>	Flip-Back Desk Arm, Swingaway detachable footrest	16"	16"*	18" or 20"	17 ¼"	29"	35½"	24"	12 ¼"	34 lb	300 lb
<b>3F012230</b>	Flip-Back Desk Arm, Elevating legrest	16"	16"*	18" or 20"	17 ¼"	29"	35½"	24"	12 ¼"	34 lb	300 lb
<b>3F012320</b>	Flip-Back Desk Arm, Swingaway detachable footrest	20"	16"*	18" or 20"	17 ¼"	29"	35½"	28"	12 ¼"	34 lb	300 lb
<b>3F012330</b>	Flip-Back Desk Arm, Elevating legrest	20"	16"*	18" or 20"	17 ¼"	29"	35½"	28"	12 ¼"	34 lb	300 lb
<b>3F012340</b>	Flip-back Full Arm, Swingaway Footrest	20"	16"*	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb
<b>3F013120</b>	Flip-Back Adjustable Height Desk Arm, Swingaway Footrest	18"	16"*	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
<b>3F013220</b>	Flip-Back Adjustable Height Desk Arm, Swingaway Footrest	16"	16"*	18" or 20"	17 ¼"	29"	35 ½"	24"	12 ¼"	34 lb	300 lb
<b>3F013320</b>	Flip-Back Adjustable Height Desk Arm, Swingaway Footrest	20"	16"*	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb

\* + 2" extension

<b>Traveler® L3 Plus – Common Accessories &amp; Parts</b>		<b>HCPCS Code</b>
<b>90763227</b>	Complete Assembly adjustable height desk arm, 1pr	E0973
<b>300-505</b>	Positioning belt 27-47" L (auto-style) push button	E0978
<b>90763062A</b>	Wheel lock lever extension, pair	E0961
<b>90763064A</b>	Heel loop, pair	E0951
<b>90763230A</b>	Anti-tippers adjustable, pair	E0971
<b>907633CSE</b>	Seat rail plug (dust cap) for seat extension, two pairs per bag	
<b>90763SE16</b>	2" seat extension kit for 16" W	E0981
<b>90763SE18</b>	2" seat extension kit for 18" W	E0981
<b>90763SE20</b>	2" seat extension kit for 20" W	E0981
<b>B10320F</b>	O2-IV holder	E2208
<b>B10325</b>	O2 cylinder holder	E2208
<b>EJO2BAG07</b>	O2 holder bag	E2208