

MEDICAL SEATING



MEDICAL SEATING

Hausted® Exam Room Furnishings line offers a variety of medical stools to meet the needs of an active facility.

Manufactured with the highest quality materials to ensure a long life, Hausted stools offer the most wanted features in an exam stool.

- Seat diameter comfortably fits to most users
- Pneumatic assist height adjustment for ease of use
- Choice of 3 base styles to complement facility décor
- Choice of height adjustment lever options
- Wide diameter base for greater stability
- Backrest and foot ring options
- Hard wheel casters
- 24 color choices



9702 Series

9702-xx physician stool with D-handle release, brushed aluminum base and toe caps



9501 Series

9501-xx physician stool with single lever release and black nylon composite base. Our most popular Physician Exam Stool — long lasting quality at an economical price



9903 Series

9903-xx physician stool with chrome 360° hand release and polished aluminum base

Specifications:

	9501 SERIES		9702 SERIES		9903 SERIES	
	inch	cm	inch	cm	inch	cm
Seat Diameter:	16.00	40.64	16.00	40.64	16.00	40.64
Base Diameter: <i>with casters</i>	24.00	60.96	24.00	60.96	24.00	60.96
Seat Height: <i>standard - low</i>	18.50	46.99	17.50	44.45	18.00 ¹	45.72 ¹
Seat Height: <i>standard - high</i>	24.00	60.96	23.00	58.42	23.00 ¹	58.42 ¹
-10BFR- Seat Height: <i>lab/counter - low</i>	24.50	62.23	23.50	59.69	22.00	55.88
-10BFR- Seat Height: <i>lab/counter - high</i>	35.00	88.90	34.00	86.36	32.00	81.28
-BFR- Overall Height: <i>with backrest - low</i>	33.00	83.82	32.50	82.55	—	—
-BFR- Overall Height: <i>with backrest - high</i>	39.00	99.06	37.50	95.25	—	—
-FR- Footring: <i>hard wheel</i>	18.00	45.72	18.00	45.72	18.00	45.72
Caster Diameter:	2.00	5.08	2.00	5.08	2.00	5.08
Weight Capacity: <i>EVENLY DISTRIBUTED</i>	lb	kg	lb	kg	lb	kg
	350	158.76	350	158.76	350	158.76
Shipping Weight:	20	9.07	22	9.98	23	10.43

Features:

Height Release: <i>hardware</i>	single lever	D-handle	360° chrome hand
Height Adjustment: <i>method</i>	pneumatic	pneumatic	pneumatic
Base Material:	black nylon composite	brushed aluminum with nylon toe caps	polished aluminum

Options:

	item number	item number	item number
None:	9501-xx	9702-xx	9903-xx
with Backrest:	9501-B-xx	9702-B-xx	—
with Footring:	9501-FR-xx	9702-FR-xx	9903-FR-xx
with Backrest & Footring:	9501-BFR-xx	9702-BFR-xx	—
Lab Height <i>with</i> Backrest & Footring:	9501-10BFR-xx	9702-10BFR-xx	—
Lab Height <i>with</i> Footring:	—	—	9903-10FR-xx

— Denotes Feature Not Available / Does Not Apply

¹ Overall Height

Standard Color Options (-xx)

Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Swatch Card is available from your representative



WARNING: This product can expose you to chemicals including Formaldehyde (gas) which is known to the State of California to cause cancer. For more information go to: www.P65Warnings.ca.gov/furniture

Color variations may occur due to monitor or printer settings.



© 2021, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com

graham field

GF1800080 Rev C21