



## CANES

# TRI-STEP CANE


The Lumex Tri-Step Cane is the ultimate cane solution designed to keep you on the go. The base provides three point contact with the ground which allows the Tri-Step to stand upright on its own and provides added stability. The contoured grip is comfortable for the hands. With its folding design and strap to secure, the cane is very convenient for storage. And to make it more eye appealing, the Tri-Step is available in 4 attractive colors.

### Features

- Three point contact base allows the cane to stand upright on its own
- Sides of base can be used to step on to lift cane if cane has fallen on the ground – no more bending over to pick up a fallen cane!
- Comfortable, contoured grip
- Lightweight aluminum, folding design makes it easy for storage
- Easily adjusts in height from 33" – 37"
- Designed to support up to 300 lb
- Product weight: 1.1 lb
- Cane tip can be purchased separately to fit most standard canes with 3/4" (19mm) tubing
- Available in Black, Blue, Green and Pink
- Designed to meet the requirements of HCPCS code: E0100
- Three (3) year limited warranty

### Replacement Parts

		UOM
6000T	Replacement Tip	1 ea



		minimum Height		maximum Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	
6000BK	Black	33	83.82	37	93.98	1.1	0.5	300	136.1	4/cs
6000B	Blue	33	83.82	37	93.98	1.1	0.5	300	136.1	4/cs
6000G	Green	33	83.82	37	93.98	1.1	0.5	300	136.1	4/cs
6000P	Pink	33	83.82	37	93.98	1.1	0.5	300	136.1	4/cs



GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company



© 2017, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.



+1.770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)

Color variations may occur due to monitor or printer settings.