



## HEALTHCARE SEATING

# SWING AWAY FOLDING TABLE

### Features

- Versatile swing-away table provides convenience to patients
- 360 degree swivel table allows items to be conveniently placed in front of the patient or to their side
- Table swings out of the way for easy patient entry or exit from recliner
- Scratch resistant table top is durable and easy to clean

★ See Chart below for Links, Item Numbers, Specifications, Options and Warranty



### Specifications

Item	Description	Meets HCPCS	Width		Length		Thickness		Height above seat		Weight Capacity <small>EVENLY DISTRIBUTED</small>		compatible recliner series						Made in USA	Warranty in Months	
			inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	FR566G	FR566DG	FR577RG	FR587W	FR587WD	FR597			FR588W
SAFT10													■	■	■	■	■	—	—	■	12
SAFT20	Swing Away Folding Table	E0184	10.00	25.40	19.00	48.26	0.75	1.905	8.50	21.59	25	1.79	—	—	—	—	—	■	—	■	12
SAFT30													—	—	—	—	—	■	—	■	12

■ Applies — Does Not Apply



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved. Lumex is a registered trademark of GF Health Products, Inc.