



LABORATORY

TEST TUBE BRUSHES

Features

• Firm white nylon bristles with radial tufted end on sturdy, twisted, galvanized wire handle

Sized to clean test tubes with a minimum of effort



	overall Length		Width		Length		
	inch	cm	inch	cm	inch	cm	UOM
3403	8	20.3	0.5	1.3	3	7.6	12 ea
3404	9	22.9	0.75	1.9	3.5	8.9	12 ea



GF Health Products, Inc. is an intertek ISO 13485:2016 / MDSAP Certified Company









