

OXIMETERS

FINGER PULSE OXIMETER

The JB02020 Finger Pulse Oximeter provides simple, fast, and accurate oxygen saturation (SpO₂) and pulse rate measurements. Its lightweight and compact design makes it an ideal solution for home, clinical and EMS settings.

Features

- Lightweight, pocket-size solution
- Provides fast and accurate SpO₂, pulse rate and Perfusion Index measurements
- Large LED digital display that is easy to read
- Simple one button operation
- Automatic power off function to conserve batteries
- Accommodates widest range of finger sizes from pediatric to adult

Specifications:

Display: Oxygen saturation (SpO₂), Pulse rate and Perfusion Index

Display Range: SpO₂: 35% to 100%

Pulse Rate: 30 to 240 bpm

Accuracy (% SpO₂): 80% - 99%: ±2% 70% - 100%: ±2%

Pulse Rate Accuracy: ±2%

Perfusion Index: 0-20%

Alarms:

SpO₂ low limit: 90%

Pulse Rate:

High limit: 120 bpm

Low limit: 50 bpm

Low Battery indicator

Unit Weight: .1 lb (52g)



	overall Length		overall Width		overall Height		UOM
	inch	mm	inch	mm	inch	mm	
JB02020	2.2	56	1.3	34	1.2	30	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
John Bunn is a registered trademark of GF Health Products, Inc.

graham
field