

## DIAGNOSTIC INSTRUMENTS

# TUNING FORKS

### Features

- Aluminum alloy tuning forks: corrosion and tarnishing resistant
- To assess a patient's hearing
- Tines are 1.75" (4.4 cm) long and 3/8" (0.95 cm) wide
- Long lasting single frequency pure tone



		overall Length		fixed weight Diameter		UOM
		inch	cm	inch	mm	
1314	C128 fixed weight	7.75	19.7	0.9	22.5	1 ea
1315	C256 fixed weight	6.5	16.5	0.7	18	1 ea
1316	C512 without weight	6.25	15.9	—	—	1 ea
1317	C1024 without weight	5.3	13.5	—	—	1 ea
1318	C2048 without weight	4.5	11.4	—	—	1 ea
1313	Set of all 5 (with cushioned vinyl zippered case)					1 set



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2012, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF1200088RevC22