



## ROLLATORS

# WALKABOUT BASIC FOUR-WHEEL ROLLATOR

The Lumex Walkabout Basic Four-Wheel Rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 13.5 lbs and supports up to a 250 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.

### Features

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable padded backbar with easy release pull knobs
- Adjustable handle height to accommodate broad range of user heights
- Pouch included for storage of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4900B	Blue	16.5	41.91	24.2	61.47	24.5	62.23	13	33.02	12.5	31.75	21	53.34
RJ4900R	Red	16.5	41.91	24.2	61.47	24.5	62.23	13	33.02	12.5	31.75	21	53.34

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4900B	Blue	6	15.24	32.5	82.55	37.5	95.25	13.5	6.12	250	113.40	1 ea
RJ4900R	Red	6	15.24	32.5	82.55	37.5	95.25	13.5	6.12	250	113.40	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2014, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

Rollators are not intended to be propelled while seated.  
Do NOT use this rollator as a wheelchair or transport device.

graham  
field

GF1400046Rev/D22