

MOBILITY ACCESSORIES

GLIDE BRAKES

Features

- For use in conjunction with swivel wheels or fixed wheels, glide brakes provide a safe, gradual braking action for rear walker legs
- Glide brakes attach to the rear of walker legs
- By pressing down on the handgrip, glide brakes allow rear walker tips to gradually touch the surface preventing forward motion
- Designed to meet the requirements of HCPCS code: E0159
- Limited lifetime warranty

Glide Brakes are recommended for use on Lumex Walkers. (We cannot guarantee the performance of these accessories on other industry products.)



		UOM
603910A	Glide Brakes	1 pr



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF2100094RevC22