

DAILY LIVING AIDS

RINSETTE

Features

- Molded styrene tray for washing hair, dressing head wounds, or irrigating ears without moving bedridden patients
- Prevents wet bed linen as excess fluids flow directly into bedside container
- Hang-up slot for convenient storage
- Cold sterilizable
- Durable and not made with natural rubber latex



		Color	UOM
5356	Heavy-Duty Rinsette	Green	1 ea
5357-1	Rinsette	Gray	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Grafcó is a registered trademark of GF Health Products, Inc.

graham
field

GF1100049RevB22