

THERAPY & REHABILITATION

MOIST HEAT THERAPY PACKS

Features

- Moist heat pack is a simple and effective method of applying moist heat to a local area
- Gives at least 30 minutes of therapeutic heat
- Use with either a moist pack cover or a heavy towel



		overall Width		overall Length		UOM
		inch	cm	inch	cm	
2088	Neck Pack	9	22.86	24	60.96	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Grafcō is a registered trademark of GF Health Products, Inc.

graham
field

GF1100007RevB22