

## SHOWER TRANSPORT CHAIRS

# PVC SHOWER BED

### Features

- Shower bed comes with closed cell water resistant foam pad
- Durable cover is easy to clean and will not promote bacteria growth
- Beveled pillow
- Three position elevating headrest
- Side safety railings are secured with user-friendly lock pins (not friction)
- Drop side railings allow barrier free transferring
- Full length drain pan catches falling water, and attached hose directs the residue to any drain, 14 drain holes
- Available in 3 anti-bacterial mesh or anti-bacterial solid vinyl colors
- Shower bed has six 4" heavy duty casters
- Center casters provide easy maneuverability and helps to eliminate sagging
- Constructed of the highest quality health care grade PVC
- Dual reinforcement into the caster fittings adds additional support
- Quick drying mesh support
- Threaded stem casters for additional safety
- Ergonomic curved rails provide support without sharp edges and reclining chairs have anti-slip handgrips for added safety
- One (1) year limited warranty on PVC frame; six (6) month limited warranty on all other materials



### Standard Mesh and Solid Vinyl Colors



Royal Blue



Forest Green



Maive

		overall Height		overall Width		overall Length		Height with safety rails		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
8005		30.5	77.47	28	71.12	79	200.66	40	101.60	450	204.12	1 ea
9005	Bariatric	31.5	80.01	32	81.28	79	200.66	40	101.60	600	272.16	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF0700130RevC22