

SHOWER TRANSPORT CHAIRS

PVC SHOWER COMMODE CHAIR

Features

- Constructed of the highest quality grade PVC
- Deluxe elongated heavy-duty open front seat
- Highest quality casters are made to last and roll easily over surfaces
- 89180 comes standard with 3" twin wheel casters; 2 locking and 2 non-locking
- 89220 comes standard with 3" single wheel casters; 2 locking and 2 non-locking
- Quick drying, removable back support is fluid resistant, anti-bacterial and anti-fungal
- One (1) year limited warranty on PVC frame; six (6) month limited warranty on all other materials



Highest quality casters are made to last



		overall Height		overall Width		overall Depth		Seat Height		Seat Depth		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
89180	3" twin wheel casters	41.3	104.90	21.5	54.61	21.5	54.61	19.5	49.53	18.3	46.48	18	45.72	300	136.08	1 ea
89220	3" single wheel casters	41.3	104.90	25.5	64.77	21.5	54.61	19.5	49.53	18.3	46.48	22	55.88	375	170.10	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2013, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-1300071RevB22